# T7 RNA Polymerase

**Code No. 2540A**      **Size:**      **5,000 U**  
**Conc.:**      **50 U/ $\mu$ l**

## Supplied Reagents:

**10X T7 RNA Polymerase Buffer**      **1 ml**  
**50 mM DTT**      **1 ml**

## Description :

T7 RNA Polymerase is an enzyme coded by the DNA of phage T7. This enzyme requires double-stranded DNA that contains the T7 promoter sequence as a template and NTPs as substrates, and synthesizes the single strand RNA complementary to the template downstream of the promoter. It has very high specificity for T7 promoter, with almost no recognition of promoters from other phages.

## Storage Buffer :

20 mM	Potassium Phosphate, pH 7.9
100 mM	NaCl
0.1 mM	EDTA
1 mM	DTT
50%	Glycerol

**Storage:**      -20°C

## Source :

*Escherichia coli* carrying a plasmid containing the gene for phage T7 RNA polymerase.

## Properties :

- Molecular Weight : approx. 98,000
- Subunit : Single polypeptide
- Optimum pH : pH 8.0
- Cofactor : Mg<sup>2+</sup> (Optimum Conc. : 8 mM)  
reductant (5 mM DTT)
- Activator : 2 mM Spermidine  
(4 mM Spermidine causes coprecipitation with DNA)
- Inhibitor : NEM (10<sup>-4</sup> M)  
PCMB (5 x 10<sup>-8</sup> M)

## Unit definition :

One unit is the amount of the enzyme that incorporates 1 nmol of [<sup>3</sup>H] GMP into the acid-insoluble fraction in 1 hour at 37°C and pH 8.0.

## Reaction mixture for unit definition :

40 mM	Tris-HCl, pH 8.0
8 mM	MgCl <sub>2</sub>
2 mM	spermidine
5 mM	DTT
0.4 mM	ATP, UTP, and CTP
0.4 mM	[ <sup>3</sup> H] GTP
1 $\mu$ g/50 $\mu$ l	pT7-2 DNA

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications :

1. Synthesis of the single strand RNA.
2. Synthesis of highly labeled RNA probes.
3. Synthesis of precursors of siRNA.
4. Synthesis of precursors for RNA splicing reactions.
5. Synthesis of capped mRNA when a cap analog is used as a primer.

**Note :** For efficient transcription of a specific region, the end of the template DNA should be blunt or should protrude in the 5'-direction downstream of the region. Spermidine contained in the buffer forms the complex with nucleic acids and may be precipitated as indissoluble materials. It is desirable to add the template DNA lastly by all means.

## Composition of Supplied Reagents (Store at -20°C) :

1. 10X T7 RNA Polymerase Buffer	
400 mM	Tris-HCl, pH 8.0
80 mM	MgCl <sub>2</sub>
20 mM	spermidine

2. 50 mM DTT

## Applications example :

10X T7 RNA Polymerase Buffer	2 $\mu$ l
50 mM DTT	2 $\mu$ l
ATP, CTP, GTP, UTP	each 2 mM
RNase Inhibitor	20 U
Template DNA	20 ng - 1 $\mu$ g
T7 RNA Polymerase	50 U
Sterile purified water	up to 20 $\mu$ l

↓  
Incubate at 42°C for 1 - 2 hr

## Related Products :

- Recombinant Ribonuclease Inhibitor (Cat. #2313A/B)
- ATP (Cat. #4041)
- GTP (Cat. #4042)
- CTP (Cat. #4043)
- UTP (Cat. #4044)
- *in vitro* Transcription T7 Kit (for siRNA Synthesis) (Cat. #6140)

## References :

- 1) Davanloo P, Rosenberg A H, Dunn J J, and Studier F W. *Proc Natl Acad Sci USA*. (1984) **81**: 2035-2039.
- 2) Chamberlin M, McGrath J, and Waskell L. *Nature*. (1970) **228**: 227-231.
- 3) Chamberlin M and Ring J. *J Biol Chem*. (1973) **248**: 2235-2244.
- 4) ET Schenborn and RC Mierendorf Jr. *Nucleic Acids Research*. (1985) **13**: 6223-6236.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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*Escherichia coli* carrying a plasmid containing the gene for phage T7 RNA polymerase

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