Human Cell-Free Protein Expression Master Mix for mRNA

Code No. 3288

Size: 180 μ l
(for 10 reactions)

Description:

Human Cell-Free Protein Expression Master Mix for mRNA is a reagent for cell-free protein synthesis using mRNA. This reagent contains all the components necessary for protein translation, including ribosomes, translation initiation/elongation factors, and tRNAs, derived from human cells. Protein synthesis is performed simply in a single tube by adding protein-encoding capped mRNA^{*1} to the reagent, followed by a few hours of incubation at 32°C. Compared to conventional cell-based experiments, this product provides much faster solutions not only for the production of target proteins but also for the evaluation of mRNA constructs. Therefore, it is suitable for the screening of mRNA constructs designed for applications such as vaccine development.

^{*1} It is important that mRNA has a cap structure for protein expression. Generally, co-transcriptional capping using cap analogs or enzymatic capping using Vaccinia or Faustovirus capping enzymes and 2'-O-methyltransferase are used for mRNA capping.

Storage: -80°C

* Store this product immediately at -80°C upon receipt. Leaving this product above -80°C may lead to a decrease in performance.

* If not used all at once, it is recommended to aliquot and store the reagent during the initial thawing to prevent performance degradation due to repeated freeze-thaw cycles.

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

1. Screening of designed mRNA constructs
2. Cell-free synthesis of toxic proteins

Protocols:

<Reaction mix>

Human Cell-Free Protein Expression Master Mix for mRNA	18 μ l
Capped mRNA	2 μ g ^{*2}
Nuclease-free water	x μ l
Total	20 μ l

1. Thaw the reagent on ice, mix well by gentle pipetting, and dispense 18 μ l into a 0.2 ml tube.
2. Pre-incubate the reagent from step 1 at 32°C for 5 minutes.
3. Add 2 μ g^{*2} of capped mRNA to the pre-incubated reagent from step 2. If necessary, add nuclease-free water to make the total reaction volume 20 μ l.
4. Mix the reaction mixture well by gentle pipetting and incubate at 27 - 32°C for 1 -5 hours.^{*3}

^{*2} Initially, 2 μ g of capped mRNA is recommended. If the production of the target protein is low, try adjusting the amount added.

^{*3} The optimum reaction temperature and time differ depending on the target protein and the mRNA construct. Initially, a reaction condition of 32°C for 3 hours is recommended.

Related Products:

Human Cell-Free Protein Expression System (Cat. #3281)
Cloning Kit for mRNA Template (BspQ I) (Cat. #6133)
Takara IVTpro™ T7 mRNA Synthesis Kit (low dsRNA) (Cat. #6134)
Takara IVTpro™ T7 mRNA Synthesis Kit (Cat. #6144)
Faustovirus Capping Enzyme (S17) (Cat. #2480A/B)
mRNA Cap 2'-O-Methyltransferase (Cat. #2470A/B)

IVTpro is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Human Cell-Free Protein Expression Master Mix for mRNA

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