

DNA MW Standard Marker  
**100 bp DNA Ladder**

**Code No. 3407A**

**Size: 500 µl**  
**(for 100 lanes)**

**Shipping at -20°C**

**Stored at -20°C**

**Supplied : 6X Loading Buffer 1 ml**  
**Reagent**

**Lot No.**

**Concentration : ng / µl**

**Volume : µl**

**Form :** 10 mM Tris-HCl, pH7.4  
1 mM EDTA

**Description :** This marker consists of 10 fragments between 100 and 1000 bp in multiples of 100 bp and an additional fragment at 1500 bp. The 500 bp band containing triple the mass of the other fragments, serves as a visible reference indicator; all other fragments appear with equal intensity on the gel. They are all double-stranded DNA fragments.

Loading volume is 5 µl per lane. 5 µl of this product contains about 50 ng of each band, except the 500 bp fragment which contains about 150 ng.

<u>Fragment</u>	<u>Size (bp)</u>
A	1,500
B	1,000
C	900
D	800
E	700
F	600
G	500
H	400
I	300
J	200
K	100

**Usage :** Used as a DNA molecular size marker in gel electrophoresis. This product can be radioactively labeled at 5'-termini by the kinase phosphorylation reaction using MEGALABEL™ (TaKaRa Cat.# 6070).

**6X Loading Buffer (Store at RT after used.) :**

36% Glycerol  
30 mM EDTA  
0.05% Bromophenol Blue  
0.05% Xylene Cyanol

Add 6X Loading Buffer with 1/5 volume of DNA solution to apply on agarose gel electrophoresis. In case precipitates generated during the storage at room temperature, dissolve in warm bath before use.

**Application example :**

100 bp DNA Ladder	5 µl
6X Loading Buffer	1 µl



Run agarose gel electrophoresis using ex. 3% NuSieve 3 : 1 Agarose (BMA Corp.).

Perform staining with EtBr, SYBR™ Green I Nucleic Acid Gel Stain or GelStar® Nucleic Acid Stain (BMA Corp.).

**Note**

For research use only. Not for use in diagnostic or therapeutic procedures.