siRNA Ladder Marker

Code No. 3430

Size: 25 μ g
(for 25 lanes)

Conc.: 200 ng/ μ l

Supplied Reagents:

6X Loading Buffer

1 ml

Description:

This marker consists of 10 dsRNA fragments between 20 and 1,000 bp, in multiple of 10 bp between 20 and 50, in multiple of 100 bp from 100 to 500 bp, and an additional fragment at 1,000 bp. The 100 bp band is observed in higher intensity to serve as a visible reference indicator.

Fragment	Size (bp)
A	1,000
B	500
C	400
D	300
E	200
F	100
G	50
H	40
I	30
J	20

Form:

10 mM Tris-HCl, pH 8.0
1 mM EDTA

Storage: -20°C

Usage:

Used as a dsRNA molecular size marker in gel electrophoresis in preparing siRNA cocktail.

6X Loading Buffer (Store at RT after used):

36% Glycerol
30 mM EDTA
0.05% Bromophenol Blue
0.02% Xylene Cyanol

Add 6X Loading Buffer with 1/5 volume of RNA solution to apply on 15% polyacrylamide gel electrophoresis. In case precipitates generated during the storage at room temperature, dissolve in warm bath (approx. 50°C) before use.


Application example:

siRNA Ladder Marker 5 μ l
6X Loading Buffer 1 μ l

Run 15% polyacrylamide gel electrophoresis.

Perform staining with EtBr, or SYBR™ Green I Nucleic Acid Gel Stain (Cat. #5760A/5761A)*.

[Electrophoresis Result] (15% polyacrylamide gel)

Note:

1. It is recommended to store in aliquots for several uses each after opened, not to repeat the excess freeze-thaw cycles.
2. Do not store siRNA Ladder Marker after adding Loading Buffer. Bands may not be observed clearly.
3. Store at -20°C quickly after used.
4. Extra precaution should be taken during the operation in order to prevent RNase contamination. Put on clean disposable gloves, and tubes and micropipette chips should be exclusively used for RNA experiments. Do not perform experiments which uses RNase in the same area.
5. Do not use this product on denatured polyacrylamide gel. The bands cannot be separated correctly since they are observed in duplicate.
6. Depending on the sequence of dsRNA sample, the fragments may appear in a little difference lengths from this ladder marker. Slight difference in size may occur when a synthesized siRNA sample partially includes DNA.

* Not available in all geographic locations. Please check for availability in your area.

SYBR is a trademark of Life Technologies Corporation.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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