EpiScope® Unmethylated HCT116 DKO gDNA

Code No. 3521

Size: 10 µg

Conc.: 100 ng/µl

* 2 years from date of receipt under proper storage conditions.

Description :

Analysis of DNA methylation is one of the regions of epigenetics research, where more attention is recently attracted. It has been known that some cytosines at CpG sequences of mammalian genomic DNA are methylated. Some CpG islands in promoter region of a gene are especially analyzed about its DNA methylation, because methylation of the promoter region may control transcriptional activity of the gene. There are various methods of DNA methylation analysis, such as Bisulfite sequence, etc. This product is genomic DNA purified from DNMT DKO (double knock-out) human HCT116 cell line, which is genetically lacking both of methyltransferases DNMT1 (DNA methyltransferase 1) and DNMT3B (DNA methyltransferase 3B). The methylation level of the genomic DNA from DNMT DKO HCT116 is extremely lower (< 5%)*¹ than the wild type genomic DNA. Therefore, the DNA can be used as a negative (non- or hypo-methylated) control*² in various experiments of DNA methylation analysis.

- * 1 This product is genomic DNA simply purified from DNMT DKO human HCT116 cell line. Therefore, the detailed methylation status of this genomic DNA could be varied between production lots.
- * 2 Please utilize EpiScope Methylated HCT116 gDNA (Cat. #3522) as a positive control (highly-methylated DNA).

Form : 10 mM Tris-HCl, pH 8.0
1 mM EDTA

Storage : -20°C

Source and preparation :

This product is genomic DNA purified from DNMT DKO HCT116 cells.

Quality Control Data :

The Methylation levels of 25 genes on EpiScope Promoter qPCR Array (Human) (Cat. #5301)*¹ were analyzed by methylated DNA enrichment using EpiXplore™ Methylated DNA Enrichment Kit (Cat. #631962) followed by qPCR analysis using 2X TB Green® *Premix Ex Taq*™ GC (Perfect Real Time) (Cat. #RR071A/B).^{*1, 2} It was confirmed that the methylation levels of all genes were low.

*1 Not available in all geographic locations. Check for availability in your area.

*2 We have begun the process of changing the names for Takara Bio's intercalator-based real-time PCR (qPCR) products to the "TB Green series". These products can be used the same way as before, as only the names are changing. Catalog number and product performance are unaffected by this transition.

Applications :

Used as a negative control (non- or hypo-methylated) human genomic DNA for Bisulfite sequence, Combined Bisulfite Restriction Analysis (COBRA), Methylated CpG Island Recovery Assay, and other DNA methylation analysis methods.

TB Green and EpiScope are registered trademarks of Takara Bio Inc.
Premix Ex Taq is a trademark of Takara Bio Inc.
EpiXplore is a trademark of Takara Bio USA, Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

EpiScope® Unmethylated HCT116 DKO gDNA

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