pAUR135 DNA

Code No. 3604

Size: 20 μ g

Conc.: 0.5 μ g/ μ l

* 2 years from date of receipt under proper storage conditions.

Description:

pAUR135 is a shuttle vector including DNA sequence to remove the vector region containing selectable marker *AUR1-C* from transformants of *Saccharomyces cerevisiae* by this vector. Selection of marker-removed clones is easy and efficient, and so pAUR135 vector is available repeatedly for transformation. The vector contains *AUR1-C* gene that confers Aureobasidin A (AbA)-resistance on yeast cells and GIN11M86 DNA under the control of GAL10 promoter. Overexpression of GIN11M86 causes growth inhibition of host cells.

Most of AbA-resistant transformants by pAUR135 vector, when shifted onto galactose medium which induce overexpression of GIN11M86 controlled by GAL10 promoter, show lethal-phenotype. But, only pAUR135 vector-removed clones as a result of low-efficient homologous recombination can grow on galactose medium preferentially. These galactose-growing clones consist of expected-phenotype ones and reverted ones, which are AbA-sensitive. So, pAUR135 is available for retransformation of these clones.

pAUR135 is effective for introduction of mutation into yeast genome and disruption of gene function. In case of transformation of industrial yeast strains, usage of pAUR135 decreases troubles with unnecessary DNA sequence.

Form: 10 mM Tris-HCl, pH 8.0
1 mM EDTA

Storage: -20°C

Preparation: cccDNA was purified by CsCl-EtBr ultracentrifugation.

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Vector of pAUR135:


AUR1-C: AbA resistant gene in *S. cerevisiae*

GAL10p-GIN11M86:

DNA sequence related to galactose-inducible growth inhibition in *S. cerevisiae*

Amp^r: Selective marker in *E. coli*

ColE1 ori: Replication origin in *E. coli*

Cloning site of pAUR135:

2245 2257 5658 5679
5'-GCATGCCTGCAGG-3' 5'-CCCGGGTACCGAGCTCGAATTC-3'
Sph I Sse8387 I Sma I Kpn I Sac I EcoR I

Protocol (disruption of gene function):

- (1) Insertion of DNA fragment for disruption into pAUR135 vector
 1. Prepare DNA fragment, whose region has a unique restriction enzyme site which pAUR135 does not, containing a longer part of the objective gene by PCR etc.
 2. Design to disrupt the targetted gene after removing pAUR135 vector with introduction of stop codons or deletion into the DNA fragment.
 3. Insert the DNA fragment to cloning site of pAUR135 and purify the plasmid vector.
 4. Linearize the plasmid vector by cutting unique site on insertion fragment.

(2) Transformation of *S. cerevisiae*

1. Transform yeast cells by acetate-lithium method.
2. Culture in 2 - 5 ml of YPD medium for more than 6 hours and then spread onto YPD selective medium plate containing AbA.
3. Confirm the disruption of targetted gene function of AbA-resistant transformants by examination of phenotype or Southern hybridization.

(3) Selection of marker (pAUR135 vector)-removed clones

1. Streak 10 - 50 clones of disruptants on YPGal plate and then isolate single colonies which can grow on galactose medium. (YPGal; 2% galactose, 2% polypeptone, 1% yeast extract)
2. Examine 10 - 100 colonies growing on YPGal plate in the AbA sensitivity and the disruption of targetted gene function.

References:

- 1) Hashida-Okado T, Ogawa A, Kato I, and Takesako K. *FEBS Letters*. (1998) **425**: 117-122.
- 2) Hashida-Okado T, Ogawa A, Endo M, Takesako K, and Kato I. *Mol Gen Genet*. (1996) **251**: 236-244.
- 3) Akada R, Yamamoto J, and Yamashita I. *Mol Gen Genet*. (1997) **254**: 267-274.
- 4) Kawahata M, Amari S, Nishizawa Y, and Akada R. *Yeast*. (1999) **15**: 1-10.

Attention:

1. This marker removing system using pAUR135 is dependent on overexpression of GIN11M86 controlled by galactose-inducible GAL10 promoter. Consequently, the galactose assimilation of host strains is required for efficient marker removal. Before using the strain for transformation by pAUR135, examine it in galactose assimilation.
2. Within marker-removed clones, ratio of expected-phenotype ones is extremely influenced by the position of introduction of mutation or deletion.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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