# dNTP Mixture

**Code No. 4030**      **Size: each 3.2  $\mu$  mol**  
**Conc.: each 2.5 mM**

\* 2 years from date of receipt under proper storage conditions.

## Description:

This product is a mixture of dATP, dCTP, dGTP, and dTTP (each 2.5 mM).

**Form:**            Solution (sodium salt), pH 7 - 9

**Storage:**        -20°C

**Purity:**            $\geq$ 98%

## Performance:

- Good performance was confirmed by PCR amplification of 500 bp and 20 kb fragments from  $\lambda$  DNA template.
- Good performance was confirmed by RT-PCR amplification of a 4.4 kb fragment of the human TFR gene from HL60 cell total RNA.

## Applications:

- 1) Ready-to-use : can be used in PCR reactions without dilution
- 2) For PCR, use :
  - 4  $\mu$ l in a 50  $\mu$ l PCR with *TaKaRa Taq*<sup>™</sup> or *TaKaRa Ex Taq*<sup>®</sup> (sufficient dNTP mix for 320 reactions)
  - 8  $\mu$ l in a 50  $\mu$ l PCR with *TaKaRa LA Taq*<sup>®</sup> (sufficient dNTP mix for 160 reactions)
- 3) Use as a substrate for other DNA amplification and cDNA synthesis reactions.

*TaKaRa Ex Taq* and *TaKaRa LA Taq* are registered trademarks of Takara Bio Inc.

*TaKaRa Taq* is a trademark of Takara Bio Inc.

### Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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