Chip for DNA Recovery TaKaRa RECOCHIP

Code No. 9039

Shipping at RT
Stored at RT

Contents :

RECOCHIP	100 pcs.
2.0 ml tube	100 pcs.
Gel cutter	3 pcs.

Lot No. (Refer to the package label)

Description :

The TaKaRa RECOCHIP is a disposable chip designed to recover DNA fragments from agarose gels without additional reagents.

Note : A microcentrifuge for 1.5 ml or 2.0 ml tubes is required to use the TaKaRa RECOCHIP.

Specification : (Refer to the Fig.1)

1. Materials

RECOCHIP :

Backing: polyester

Non-woven fabric: polyester

Dialysis membrane : cellulose

(Nominal fractionating MW : 12,000-14,000)

2.0 ml tube : polypropylene

Gel cutter : polyester

2. Dimensions

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
2.0 ml tube : 39 mm length, 9 mm diameter *

(* : inside of tube)

Gel cutter : 35 (L) × 10 (W) mm

Fig. 1 RECOCHIP

Protocol : (Refer to Fig. 2)


- After electrophoresis, transfer the stained agarose gel onto UV transilluminator (or use hand illuminator) to confirm the DNA band to recover. Push the gel cutter into the agarose gel on the anode (+) side near the DNA band to create a gap.
- Dip the RECOCHIP in electrophoresis buffer to fill up the recovery space (the space between the black surface (non-woven fabric) and the dialysis membrane) of the RECOCHIP with buffer.
- Push the RECOCHIP into the gap. Make sure that the black surface (non-woven fabric) of the RECOCHIP faces toward the DNA band and the white surface (dialysis membrane) faces toward the anode. * 1
- Return the agarose gel to the electrophoresis chamber and perform electrophoresis again. * 2
- Stop electrophoresis and remove the RECOCHIP from the agarose gel. Put the RECOCHIP into the provided 2.0 ml tube and centrifuge at 5,000 rpm for 5 sec.
- 30-50 μ l of DNA solution will be recovered in the 2.0 ml tube. Perform ethanol precipitation to collect the DNA if necessary.

* 1 : Make sure the direction of the RECOCHIP is correct. DNA cannot be recovered if the direction of RECOCHIP is reversed.

* 2 : Loading time varies depending on the size of the DNA fragment, agarose gel concentration, and voltage. Please refer to the table below to determine the appropriate loading time. Takara recommends confirming the mobility and transfer of the DNA band by checking on a UV illuminator again after recovery.

DNA fragment size (Amount)	Gel conc.	Buffer	Voltage	Loading time
1.88 kb (39 ng)	1.2%	TAE	100 V	5 min.
1.88 kb (39 ng)	0.7%	TAE	100 V	3 min.
1.88 kb (39 ng)	0.7%	TAE	50 V	6 min.
1.49 kb (30 ng)	1.0%	TAE	100 V	2 min.
0.33 kb (7 ng)	2.5%	TAE	100 V	2 min.

Fig. 2 Directions for RECOCHIP use

Notice for :

- Do not autoclave the RECOCHIP.
- Do not treat the RECOCHIP with solvents, including chloroform.
- Use tweezers or wear gloves when handling the RECOCHIP. Do not touch the RECOCHIP directly with bare hands.

Note

This product is intended to be used for research purpose only. They are not to be used for drug or diagnostic purposes, nor are they intended for human use. They shall not to be used products as food, cosmetics, or utensils, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please call at +81 77 543 7247 or contact from our website at www.takara-bio.com.

Manufactured by Toyo Roshi Kaisha, Ltd.

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