For Research Use

TaKaRa

TaKaRa BioMasher Standard (Non-sterile)
TaKaRa BioMasher Standard (Sterile)

Product Manual

Table of Contents

l.	Description	3
II.	Components	3
III.	Storage	3
IV.	Materials	3
V.	Protocol	3
VI.	Experimental Examples	4
VII	Related Products	5



I. Description

TaKaRa BioMasher Standard is a disposable microtube homogenizer designed to efficiently crush small amounts of biological samples for nucleic acid and/or protein extraction. This product is a set of a 1.5-ml microtubes and micro stir bars. The inner wall of the tube and the tip of the stir bar are textured allowing effective disruption of the sample. In addition, a lid in the middle of the micro stir bar can be prevented scattering of samples and reagents during crushing. After the sample has been homogenized, the tube can be centrifuged to minimize sample loss.

The TaKaRa BioMasher Standard is available either sterilized (Cat. #9791A) or nonsterilized (Cat.# 9790A).

II. Components

- 1. TaKaRa BioMasher Standard (micro stir bar and microtube)
- 2. Stir bar grip (PESTLE GRIP)



50

III. Storage

Room temperature

* Keep out of sun and UV light.

IV. Materials

TaKaRa BioMasher Standard

Polypropylene microtubes Autoclavable
Polyacetal stir bar Not autoclavable*
Silicon rubber stir bar grip Autoclavable

* The stir bar can not be autoclaved. If you require sterilization, please use Cat. #9791A.

V. Protocol

NOTE: Wear a mask, gloves, lab coat, and safety goggles. If necessary, perform in a safety cabinet.

- 1. Insert the micro stir bar into the end of the PESTLE GRIP.
- 2. Add the sample to the microtube. Use less than 100 mg of sample. If necessary, add the extraction reagent (less than 250 μ l).
- 3. Insert the micro stir bar into the microtube. With the PESTLE GRIP, rotate and move the micro stir bar up and down, crushing the sample. If necessary, perform on ice.
- 4. After crushing, disconnect the PESTLE GRIP and discard the micro stir bar. The PESTLE GRIP can be used repeatedly; store for future use.
- 5. Close the lid of the microtube. The sample is ready for extraction.

VI. Experimental Example


[Total RNA extraction from mouse liver]

100 mg of mouse liver was crushed using the TaKaRa BioMasher Standard (Sterile). As a comparison, 100 mg frozen mouse liver was crushed with a mortar or with bead crusher. All samples were extracted total RNA using 1 ml RNAiso Plus (Cat. #9108/9109) according to the recommended protocol. Total RNA was quantified and analyzed on an 1% agarose gel.

[Results]

Processing the samples using the TaKaRa BioMasher Standard had similar extraction efficiency as the mortar and bead crushing protocols.

Crushing method	total RNA (μ g)	A_{260}/A_{280}
1. Mortar	360	1.70
2. Mortar	424	1.80
3. Bead crusher	530	1.92
4. Bead crusher	494	1.93
5. TaKaRa BioMasher Standard	411	1.89
6. TaKaRa BioMasher Standard	431	1.90



M: Wide-Range DNA Ladder (100 - 2,000 bp)

1-6: Equal amounts of RNA

TaKaRa BioMasher Standard (Non-sterile) TaKaRa BioMasher Standard (Sterile)

Cat. #9790A/9791A **TakaRa**

VII. Related Products

RNAisoPlus (Cat. #9108/9109) NucleoSpin RNA (Cat. #740955.10/.50/.250) NucleoSpin RNA XS (Cat. #740902.10/.50/.250) NucleoSpin Tissue (Cat. #740952.10/.50/.250) NucleoSpin Tissue XS (Cat. #740901.10/.50/.250)

NOTE: This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.