

For Research Use

---

**TAKARA**

***Brevibacillus* Competent Cells**

---

Product Manual

---

## Table of Contents

I.	Description.....	3
II.	Components .....	3
III.	Materials Required but not Provided .....	3
IV.	Storage .....	3
V.	Protocol.....	4
VI.	Transformation Efficiency.....	6
VII.	Genotype.....	6
VIII.	Related Products .....	6
IX.	Notice: Living Modified Organism.....	6

## ***Brevibacillus* Competent Cells**

### **I. Description**

*Brevibacillus* Competent Cells are host cells used for transformation with *Brevibacillus* expression vectors such as pNY326 DNA (Cat. #HB111), pNCMO2 DNA (Cat. #HB112), pNC-HisT DNA (Cat. #HB121), and others. *Brevibacillus* Competent Cells are prepared for transformation using the New Tris-PEG (NTP) method for chemical transformation of *B. choshinensis* strain SP3. With the NTP method, it is possible to obtain a transformation efficiency that is approximately equal to that obtained with electroporation. The NTP method also enables direct transformation of a DNA ligation solution without DNA purification steps including ethanol precipitation. For instructions on the use of the *Brevibacillus* expression system, please refer to the manual for BIC System (Cat. #HB300) or *Brevibacillus* Expression System II (Cat. #HB200) for secretory protein expression and the manual for pNI DNA/pNI-His DNA (Cat. #HB131/ HB132) for intracellular protein expression.

### **II. Components**

<i>Brevibacillus</i> Competent Cells	100 $\mu$ l x 10
MT Medium	1 ml x 10
Solution A	1 ml
Solution B	1 ml x 2

### **III. Materials Required but not Provided**

#### **1. Reagents**

- Plasmid for target gene expression
- MTNm plates
- MT Liquid Medium
- Neomycin

**Note:** See Section V-1. for MTNm plate and MT Liquid Medium components and preparation.

#### **2. Materials**

- Sterile culture tubes
- Sterile microtubes
- Microcentrifuge
- Vortex
- Incubating orbital shaker

### **IV. Storage** -80°C

**Note:** Please store at -80°C or less. Insufficient temperature control may lead to decreased transformation efficiency. Do not store in liquid nitrogen.

**V. Protocol****1. Preparation**

Prepare the following reagents and materials:

Plasmid for target gene expression  
MTNm Plates\*<sup>1</sup>  
MT Liquid Medium\*<sup>1</sup>  
Sterile culture tubes  
Sterile microtubes

**\*1 Medium Composition**

MT Liquid Medium	
Glucose* <sup>2</sup>	10.0 g/L
Phytone Peptone	10.0 g/L
35%Ehrlich Bonito Extract	5.75 g/L
Yeast extract Blue label	2.0 g/L
FeSO <sub>4</sub> · 7H <sub>2</sub> O	10 mg/L
MnSO <sub>4</sub> · 4H <sub>2</sub> O	10 mg/L
ZnSO <sub>4</sub> · 7H <sub>2</sub> O	1 mg/L
MgCl <sub>2</sub> · 6H <sub>2</sub> O	4.1 g/L
Adjust to pH 7.0 with NaOH	

\*<sup>2</sup> Sterilize glucose and glucose-free media separately. Mix after sterilization.

**MTNm Plates**

Suspend 7.5 g of agar in 500 ml of MT Liquid Medium and sterilize using an autoclave. Let stand at room temperature until it has cooled to approximately 50°C and then add neomycin solution (50 mg/ml stock solution) to a final concentration of 50 μg/ml. Mix gently and dispense into plates.

For the components of the MT Medium, the following manufacturers are recommended.

Phytone Peptone	(Becton Dickinson, Code. 211906)
35%Ehrlich Bonito Extract	(Kyokuto Pharmaceutical, Code. 551-01212-5)
Yeast extract Blue label	(Oriental Yeast Co., Ltd.)

## ***Brevibacillus* Competent Cells**

---

### **2. NTP Transformation Method**

- (1) Thaw Solution A, Solution B, and MT Medium.
- (2) Transfer only the number of tubes of *Brevibacillus* Competent Cells needed for transformation from storage, and keep on dry ice/ethanol.
- (3) Thaw the *Brevibacillus* Competent Cells quickly (approximately 30 seconds) in a 37°C water bath.
- (4) Centrifuge the cells (12,000 rpm, 30 seconds to 1 minute) to form a cell pellet and remove the supernatant with a micropipette.

Perform the following procedures at room temperature.

- (5) Mix the plasmid DNA solution (in a volume of 5  $\mu$ l or less)\*<sup>1</sup> with 50  $\mu$ l of Solution A.
- (6) Add all of the DNA solution to the *Brevibacillus* cell pellet (from step 4) and vortex to completely suspend the pellet.
- (7) Allow to stand for 5 minutes at room temperature.
- (8) Add 150  $\mu$ l of Solution B (PEG solution)\*<sup>2</sup> and vortex until the solution is uniform (5 - 10 seconds).
- (9) Centrifuge the cells (5,000 rpm, 5 minutes) and remove the supernatant.
- (10) Centrifuge briefly (5,000 rpm, 30 seconds) and remove the supernatant completely.
- (11) Add 1 ml of MT Medium and suspend completely with a micropipette.
- (12) Incubate at 37°C in an orbital shaker (120 rpm, 2 hours).
- (13) Use a sterile inoculating loop to remove a small sample from the culture. Streak on the MTNm plates and culture overnight at 37°C.
- (14) Select isolated colonies for plasmid analysis or protein expression.

\*1 When DNA ligation solution is used, mix 5  $\mu$ l of the reaction solution with Solution A. When using purified plasmids, use 10 - 100 ng.

\*2 Solution B (PEG solution) is highly viscous - use a 1,000  $\mu$ l micropipette and pipette slowly.

## **VI. Transformation Efficiency**

Transformation with 10 ng of pNY326 plasmid was performed according to the protocol and colonies that formed on a MTNm plate were selected. Transformation efficiency was  $>10^5$  transformants/ $\mu$ g pNY326 plasmid.

## **VII. Genotype**

An essential gene for spore formation is disrupted in *B. choshinensis* SP3; therefore sterilization of transformants may be performed using standard autoclave conditions. Furthermore, trace activity of an intracellular protease (*imp*) and extracellular protease (*emp*) have been disrupted to prevent degradation of expressed proteins.

**VIII. Related Products**

BIC System (Cat. #HB300)  
*Brevibacillus* Expression System II (Cat. #HB200)  
pNY326 DNA (Cat. #HB111)  
pNCMO2 DNA (Cat. #HB112)  
pNY326-BLA DNA (Cat. #HB114)  
pNC-HisT DNA (Cat. #HB121)  
pNC-HisF DNA (Cat. #HB122)  
pNC-HisE DNA (Cat. #HB123)  
pNI DNA (Cat. #HB131)  
pNI-His DNA (Cat. #HB132)

**IX. Notice: Living Modified Organism**

*Brevibacillus* Competent Cells (Cat. #HB116), BIC system (Cat. #HB300), and *Brevibacillus* Expression System II (Cat. #HB200) include a genetically "Living Modified Organism (LMO)" defined in "The Cartagena Protocol on Biosafety". The supplied *Brevibacillus* Competent Cells in these kits contain partial sequences of the 2  $\mu$ m plasmid derived from *Saccharomyces cerevisiae*.

Please follow the guidelines, laws, and regulations specific to your country and ensure safe handling, storage, transport, and disposal.

This product is developed and manufactured by Higeta Shoyu Co., Ltd. and sold by Takara Bio Inc.

**NOTE:** This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takara-bio.com](http://www.takara-bio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.