

For Research Use

---

# TaKaRa

**PrimerArray<sup>®</sup> Analysis Tool for  
Hepatic Differentiation (Human)**

---

Manual

---

## Table of Contents

I. Calculating and exporting Ct values .....	3
II. Relative quantification analysis .....	4
III. Troubleshooting .....	9

The PrimerArray Analysis Tool for Hepatic Differentiation (Human) is a software tool for analysis of data obtained using PrimerArray Hepatic Differentiation (Human) (Cat. #PH017), a primer set for real-time RT-PCR analysis of gene expression related to hepatic differentiation. The tool allows comparison between data for an unknown and control sample and performs relative quantitative analysis using Ct values exported from real-time PCR instrument software by the  $\Delta\Delta$  Ct method. Results are displayed in a graphical format.

- \* The PrimerArray Analysis Tool for Hepatic Differentiation (Human) uses a Microsoft Office Excel format file containing macros. Its performance has been validated in the following operating systems and versions of Microsoft Office Excel:
  - Windows XP operating system
  - Microsoft Office Excel 2003
  - Microsoft Office Excel 2007
- \* The PrimerArray Analysis Tool for Hepatic Differentiation (Human) is available for download from the Takara Bio website.

## I. Calculating and exporting Ct values


Set the analysis parameters using the real-time PCR instrument software, and calculate Ct values. Refer to the instruction manual of the real-time PCR analysis software for specific details of the analysis procedure.

### (1) Setting analysis parameters


The analysis parameters are automatically set in most real-time PCR analysis software. However, settings should be reviewed to ensure that those parameters are correct. If they are incorrect, the parameters will need to be re-set manually.

#### Baseline region

Set the flat region before the amplification curve begins to rise as the baseline region. If this region is not long enough, the baseline will not be properly normalized. In contrast, if this region is too long, it may cause the amplification curve to lower progressively (refer to the graphs below).




Correct Baseline



Baseline Region is Too long

#### Threshold

Set the threshold within the region of exponential PCR amplification. This is the region where the amplification curve becomes linear when vertical axis of the curve is plotted on a log scale.



Correct Threshold

- (2) Calculation of Ct value  
Most real-time PCR analysis software automatically calculates the Ct value.
- (3) Output of the data  
Output of the Ct values is generally in Microsoft Office Excel or CSV format. The output form varies depending on the analysis software used.  
\* Some real-time PCR analysis software packages do not output data from wells where sample information is not set or from wells omitted from the analysis. In this case, errors are likely during data input into the PrimerArray Analysis Tool for Hepatic Differentiation (Human). Ensure data from all wells is exported before using the analysis tool.

## II. Relative quantification analysis

Below is a protocol to perform relative quantitative analysis using the  $\Delta\Delta$  Ct method with the PrimerArray Analysis Tool for Hepatic Differentiation (Human).

- (1) Starting the PrimerArray Analysis Tool for Hepatic Differentiation (Human)  
Open the PrimerArray Analysis Tool for Hepatic Differentiation (Human) (PrimerArray Analysis Tool for Hepatic Differentiation (Human).xls) file.
- (2) Select a plate  
Choose the PrimerArray plate used for your experiment, then click the "Plate Select" button.

PlateList		
Human	Product Code	Product Name
+	PH017	PrimerArray Hepatic Differentiation (Human)

**Plate Select**

- (3) Input Control Sample Data  
After clicking "Plate Select" button, a sheet for control sample data will appear. Input Ct values in exp1 (C column), exp 2 (D column), exp 3 (E column), etc. This can generally be done by copying and pasting the Ct value output from the real-time PCR analysis software. Data for up to 10 repeated experiments can be entered.

	A	B	Control Sample										M	N
	Symbol	Well	exp1	exp2	exp3	exp4	exp5	exp6	exp7	exp8	exp9	exp10	AVG	SD_Ct
3	NANOG	A01	20.58										20.58	0.00E+00
4	GRB7	A02	23.24										23.24	0.00E+00
5	GSC	A03	25										25.00	0.00E+00
6	CYP3A7	A04	---										0.00	0.00E+00
7	ASGR1	A05	33.75										33.75	0.00E+00
8	TAT	A06	32.15										32.15	0.00E+00
9	NR1I2	A07	27.87										27.87	0.00E+00
10	CYP2A6	A08	30.69										30.69	0.00E+00
11	ABCC3	A09	32.1										32.10	0.00E+00
12	UGT1A1	A10	33.14										33.14	0.00E+00
13	ASS1	A11	22.38										22.38	0.00E+00
14	GUSB	A12	23.98										23.98	0.00E+00
15	POU5F1	B01	17.55										17.55	0.00E+00
16	IFITM1	B02	18.84										18.84	0.00E+00
17	NODAL	B03	23.65										23.65	0.00E+00
18	DLK1	B04	27.7										27.70	0.00E+00
19	KRT8	B05	32.86										32.86	0.00E+00
20	TF	B06	29.71										29.71	0.00E+00
21	NR1H3	B07	---										0.00	0.00E+00
22	CYP2C8	B08	30.1										30.10	0.00E+00

Navigation tabs: PlateSelect, PlateInfo, GeneID, TestSampleData, **ControlSampleData**, PCR\_amp\_eff

(4) Input Test Sample Data

Select the sheet "TestSampleData" for Test Sample data input. Input the data in the same way as the Control Sample. After inputting the data, click the "set sample data" button.

Clearing data

If you need to re-input data, click the "clear" button. This will delete all of the data.

Setting the Ct value cutoff

Once a Ct value cutoff is set, Ct values beyond a certain level will be excluded from analysis. The default cutoff is set at 35 cycles, and will exclude Ct values greater than 35. To change this cutoff level, change the "Ct cutoff value" .

(5) Calculation of the Normalization Factor

Click on "Set Sample Data" . The sheet "normalization\_factors" should open for calculation of the Normalization Factor. Select housekeeping gene (HKG)\*1 for normalization by checking the box in the column A, and then clicking the "NF value" button. The Normalization Factor is calculated and relative quantitative analysis will be performed automatically.

HKG	Test Sample		Control Sample		NF value
	Quantity	SD_Q	Quantity	SD_Q	
<input checked="" type="checkbox"/> GUSB	1.01E-07	0.00E+00	6.04E-08	0.00E+00	
<input checked="" type="checkbox"/> HPRT1	2.11E-08	0.00E+00	2.82E-07	0.00E+00	
<input checked="" type="checkbox"/> PGK1	1.21E-07	0.00E+00	6.70E-07	0.00E+00	
<input checked="" type="checkbox"/> ACTB	2.70E-06	0.00E+00	1.61E-05	0.00E+00	
<input checked="" type="checkbox"/> GAPDH	2.19E-06	0.00E+00	9.20E-06	0.00E+00	
<input checked="" type="checkbox"/> TBP	7.61E-09	0.00E+00	7.60E-08	0.00E+00	
<input checked="" type="checkbox"/> B2M	1.64E-06	0.00E+00	2.09E-07	0.00E+00	
<input checked="" type="checkbox"/> PPIA	1.25E-06	0.00E+00	5.74E-06	0.00E+00	
Normalization factors	Quantity	SD_Q			
NF Test	2.67E-07	0.00E+00			
NF Control	7.92E-07	0.00E+00			

\* 1 Selection of housekeeping gene:

The normalization factor is the coefficient used to normalize the template quantities used in the reaction. A housekeeping gene (HKG) whose expression level is stable among the samples is used as the index for this calculation. Care should be taken in selecting the housekeeping gene, because incorrect results can be obtained if a gene having differing expression levels among samples is used as an index. To select an appropriate housekeeping gene, confirm stable expression experimentally or use known information (biological insight, published literature, microarray analysis results, etc.).

If there is no known information suggesting an appropriate gene, use all of the control housekeeping genes as a reference. Alternatively, perform the analysis without normalization of the RNA amount (without Housekeeping Gene Normalization).

**References**

- Housekeeping Gene Primer Set (Cat. #3790/3791/3792)\*2
- geNorm manual  
[http://medgen.ugent.be/~jvdesomp/genorm/geNorm\\_manual.pdf](http://medgen.ugent.be/~jvdesomp/genorm/geNorm_manual.pdf)
- Vandesompele J, *et al.* Accurate normalization of real-time quantitative RT-PCR data by geometric averaging of multiple internal control genes. *Genome Biol.* (2002) Jun18; **3**(7): RESEARCH0034. Epub 2002 Jun 18.

\* 2 Not available in all geographic locations. Check for availability in your area.

- (6) Confirmation of the analysis results  
After the analysis, a 3D profile of the Fold Differences will appear. Select each sheet to view the additional results.

Fold Difference


The list will show the relative quantification values (fold difference) and standard deviation of the Test Sample, with the Control Sample set to 1. Standard deviation is only displayed if repeated experiments have been performed.

Fold Difference					
Well	Symbol	expression level		SD	
		Test Sample	Control Sample	Test Sample	Control Sample
A01	NANOG	0.02	1.00E+00		
A02	GRB7	0.16	1.00E+00		
A03	GSC		1.00E+00		
A04	CYP3A7				
A05	ASGR1	3.08E+04	1.00E+00		
A06	TAT	1.12E+04	1.00E+00		
A07	NR1I2	21.11	1.00E+00		
A08	CYP2A6	75.38E+02	1.00E+00		
A09	ABCC3	989.12	1.00E+00		
A10	UGT1A1	67.47E+02	1.00E+00		
A11	ASS1	20.11	1.00E+00		
A12	GUSB	4.96	1.00E+00		
B01	POU5F1	0.00	1.00E+00		
B02	IFITM1	0.26	1.00E+00		
B03	NODAL		1.00E+00		
B04	DLK1	0.11	1.00E+00		
B05	KRT8	6.63	1.00E+00		
B06	TF	2.89E+04	1.00E+00		
B07	NR1H3				
B08	CYP2C8	63.83E+02	1.00E+00		
B09	ABCG2	7.67	1.00E+00		
B10	GSTA2	1.29E+04	1.00E+00		
B11	CPS1	132.51	1.00E+00		
B12	HPRT1	0.22	1.00E+00		

Scatter plot

The left table shows a list of values and standard deviations before relative quantification with the Control Sample. The values are shown in Scatter plot in the graph on the right.

Well	Symbol	expression level		SD	
		Test Sample	Control Sample	Test Sample	Control Sample
		A01	NANOG	1.49E-02	8.06E-01
A02	GRB7	2.00E-02	1.28E-01		
A03	GSC		3.76E-02		
A04	CYP3A7	1.54E-01			
A05	ASGR1	2.69E+00	8.74E-05		
A06	TAT	2.97E+00	2.65E-04		
A07	NR1I2	1.09E-01	5.15E-03		
A08	CYP2A6	5.50E+00	7.29E-04		
A09	ABCG2	2.71E-01	1.74E-04		
A10	UGT1A1	9.00E-01	1.38E-04		
A11	ASS1	4.65E+00	7.01E-01		
A12	GUSB	3.79E-01	7.60E-02		
B01	POU5F1	1.05E-02	6.50E+00		
B02	IFITM1	1.11E-01	7.69E+00		
B03	NODAL		9.60E-02		
B04	DLX1	6.38E-04	5.79E-09		
B05	KRT8	1.08E-03	1.62E-04		
B06	TF	4.16E+01	1.44E-02		
B07	NR1H3	3.48E-01			
B08	CYP2C8	7.07E+00	1.10E-02		
B09	ABCG2	6.20E-02	8.25E-01		
B10	GSTA2	7.09E+00	5.15E-04		
B11	CPST	6.18E+00	4.67E-02		
B12	HPRT1	7.90E-02	3.96E-01		
C01	SOX2	8.79E-04	8.50E-01		
C02	PODXL	1.66E-01	2.29E+00		
C03	FOXA2	1.53E-01	2.21E-04		
C04	PROX1	7.87E-01	4.45E-01		



3D Profile

A table listing the Fold Difference of the Test Sample and gene symbols is shown, with the placement of the data corresponding to the position on the plate. The color is indicative of the degree of expression difference: red, increased expression (fold difference > 2); gray, minimal change (fold difference 0.5 - 2); blue, reduced expression (fold difference < 0.5); yellow, Ct value is not detected in both samples or one sample.


3D Profile	1	2	3	4	5	6	7	8	9	10	11	12
A	1.85E-02	1.57E-01			3.88E+04	1.12E+04	2.11E+01	7.54E+03	9.83E+02	6.15E+03	2.81E+01	4.39E+08
B	1.69E-03	2.84E-01		1.10E-01	6.63E+00			6.39E+03	7.87E+00	4.29E+04	1.33E+02	2.22E-01
C	1.83E-02	6.94E-03	6.63E+01	3.64E+03	1.75E+00			6.57E+04		1.14E+00	3.29E-01	5.36E-01
D	2.82E-01	1.30E-02	3.32E+00	2.78E+02	8.97E+04	1.84E+01	6.77E+00	5.99E+02	5.11E+04	6.93E-01	3.41E+01	4.97E-01
E	6.91E-04		1.01E+01	3.56E+04	5.24E+09	9.78E+00	5.11E+04			1.39E+01	5.46E-02	7.07E-01
F	7.48E-03	7.41E+00	8.56E+01			7.19E+02	1.70E+01	1.16E+03		6.97E-01	1.40E+00	2.97E-01
G		1.26E-01	1.57E+01		1.52E+01	7.62E+00	4.10E+02	3.53E+02			4.91E+02	2.33E+01
H		5.73E+01			8.30E+04	6.19E+00	3.25E+04	8.72E-02	8.31E+03	1.98E+01	5.17E+00	6.97E-01

Symbol	1	2	3	4	5	6	7	8	9	10	11	12
A	NANOG	GRB7	GSC	CYP3A7	ASGR1	TAT	NR1I2	CYP2A6	ABCG2	UGT1A1	ASS1	GUSB
B	POU5F1	IFITM1	NODAL	DLX1	KRT8	TF	NR1H3	CYP2C8	ABCG2	GSTA2	CPST	HPRT1
C	SOX2	PODXL	FOXA2	PROX1	KRT8	TDO2	NR1H4	CYP2E1	SLC10A1	ATP5G1	NAG5	PGK1
D	CD9	TDO2	SOX17	IFITM1	EGG	NR1A	NR1A	CYP2C9	SLC22A1	POU5F1	ACT1	ACT1
E	DMRT3B	ZFP42	CXCR4	SEPPORNA1	EGG	NR1A	NR1C1	CYP3A4	SLC22A2	PPARGC1A	KRT19	GAPDH
F	GABRB3	SOX7	GATA4	AHSG	G6PC	NR1A	PPARA	CYP7A1	SLC10B1	UCP2	KRT17	TBP
G	GAL	LAMB1	GATA6	ALB	PCSK2	AHR	CYP1A1	ABCB4	SLC10B3	ARG1	RBP1	NR2F1
H	GDF3	SRFB	AFP	APOM4	BSRP	NR1D3	CYP1A2	ABCC1	SLC20B1	ANS	HEH1	PPA

Specify a Category, then click the “Reload” button. Changes in expression of genes in that category are displayed in the table, with the color corresponding to the degree of expression difference, and as a bar graph. To display all of the data, click “All” to select all of the categories. Click “Clear” to remove all selected categories.

		All	Clear	Reload
Category				
<input type="checkbox"/>	Pluripotent/ES Cell Marker			
<input type="checkbox"/>	Extraembryonic Endoderm			
<input type="checkbox"/>	Mesendoderm			
<input type="checkbox"/>	Endoderm			
<input type="checkbox"/>	Early Hepatic			
<input checked="" type="checkbox"/>	Hepatic			
<input checked="" type="checkbox"/>	Hepatic Nuclear Receptors			
<input type="checkbox"/>	CYP Enzymes			
<input type="checkbox"/>	Hepatic Transporters			
<input type="checkbox"/>	Conjugating Enzymes			
<input type="checkbox"/>	Mitochondrial Genes			
<input type="checkbox"/>	Urea Cycle			
<input type="checkbox"/>	Cholangiocyte/Biliary			
<input checked="" type="checkbox"/>	House Keeping Genes			



Analysis is complete. When continuing the analysis with a different data set, erase the data by clicking the “clear” button on the “TestSampleData” sheet. Begin again at step (2) Select a Plate.



### III. Troubleshooting

- Security alert appears.  
PrimerArray Analysis Tool for Hepatic Differentiation (Human) includes a macro, and a security alert may appear. In this case, enable macros

Microsoft Office Excel 2007

- (1) Click “Options” on the security warning.



- (2) Select the “Enable this content”, and then click the OK button.



PrimerArray is a registered trademark of TAKARA BIO INC.

**NOTE:** This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.