

Lonza Rockland, Inc. www.lonza.com scientific.support@lonza.com Scientific Support: 800-521-0390 Document # 00525284-0911-01 Rockland, ME 04841 USA

## ProSieve<sup>™</sup> EX transfer buffer

Instructions for use

### Introduction

ProSieve<sup>™</sup> EX transfer buffer for Western blotting is a modified Towbin classic buffer. This proprietary formulation drastically increases the transfer speed without excessive heat generation affecting the proteins.

ProSieve<sup>™</sup> EX transfer buffer does not require methanol or SDS, but it is compatible with these additives to slow the higher molecular weight protein transfer.

#### Contents

Supplied as a 10X concentrate

Store at room temperature for up to 1 year from the date of manufacture.

### Protocol

- 1. Prepare a 1X transfer buffer solution by diluting the stock buffer (10X solution) with distilled water.
- 2. Soak the gel for 5 minutes in distilled water, then 5 minutes in 1X ProSieve<sup>™</sup> EX transfer buffer.
- 3. Soak the membranes and pads for at least 5 minutes in 1X ProSieve<sup>™</sup> EX transfer buffer.

Note: if necessary pre-wet transfer membrane prior to equilibration following your normal protocol (e.g. methanol for PVDF)

4. Prepare the transfer assembly as usual for your tank or semi-dry blot following the manufacturer instructions.

Suggested running conditions for tank and semidry assemblies:

Tank blot: Current transfer 375 mA

Time 10-15 min

Semidry blot:	
Constant settings	Time
375 mA	10-15 min
or	
25 V	10-15 min

#### Notes

 Buffer is compatible with all suitable membrane types used in western blotting (PVDF, Nitrocellulose, Nylon, etc...)

#### **Product Safety**

For details regarding product safety, see material safety data sheet (MSDS); call +1 (800) 638-8174 for extra copies of the MSDS. Emergency after hours, call collect +1 (303) 595-9048.

#### Ordering Information:

Catalog #	Description	Size
00200309	ProSieve™ EX transfer buffer	1L

#### **Related Products:**

PAGEr<sup>™</sup> gold tris-glycine precast gels PAGEr<sup>™</sup> minigel chamber and blot module kit ProSieve<sup>™</sup> EX running buffer ProSieve<sup>™</sup> ProTrack<sup>™</sup> dual color protein loading buffer ProSieve<sup>™</sup> blue protein staining solution ProSieve<sup>™</sup> protein markers (stained and unstained) AccuGENE<sup>™</sup> molecular biology grade water

#### For research use only. Not for use in diagnostic procedures.

# Buffer technology licensed from DGel Electrosystem, patents pending

©2011 Lonza Rockland, Inc. All rights reserved.