

β Gal-Juice/ β Gal-Juice PLUS (1 033 11/ 1 033 12)

Components included:

β-Gal-Juice/plus

50ml

Chemiluminescent Substrate for measurement of β-Galactosidase Enzyme.

Store at +4°C.

Triggering Reagent

50ml

Store at +4°C

Preparation of Cell Lysates:

For cell lysis, we suggest the follow buffers:

- **Animal cells:** 100 mM potassium phosphate (pH 7,8), 0,2% Triton X100
- **Yeast cells :** 0,5 M sodium phosphate (pH 7,1) 50 mM KCl, 5 mM MgSO₄

Standard Protocol for Detection of β-Galactosidase in Microtiter Plate Luminometer


Equilibrate at room temperature for 30 minutes before use.

Best results for chemiluminescent detection of β-Galactosidase enzyme or reporter assays can be obtained from 30 minutes to 60 minutes incubation of β-Gal Juice with β-Galactosidase enzyme and read the plate or tube within 2 to 10 minutes after triggering the reaction mixture by accelerator or enhancer.

- Add 5 to 10 µl of diluted β-Galactosidase enzyme or cell extract to microplate wells.
- Add 50 to 100 µl of β-Gal Juice and incubate for 30 to 60 minutes at room temperature.
- Add 50 to 100 µl of triggering reagent and shake it for 10 seconds.
- Place the microtiter plate in luminometer and start reading.
- Best results are obtained between 2 to 10 minutes.

β Gal-Juice/ β Gal-Juice PLUS (1 033 11/ 1 033 12)

Standard procedure:



Optional: parallel measurement in Duo- or Triple-Well System (Firefly/ Renilla/ Gaussia/ AP/ βGal)

Split the cell lysate and transfer to different tubes for measuring several enzymes in one sample.

