# PrimeSTAR® Max DNA Polymerase Ver.2

**Code No. R047S**      **Size: 625  $\mu$ l**  
**(for 25 reactions)**

**Code No. R047A**      **Size: 625  $\mu$ l x 4**  
**(for 100 reactions)**

## Description:

PrimeSTAR Max DNA Polymerase Ver.2 is a high-performance PCR amplification system that offers the highest levels of fidelity and extension speed. By adding factors that enable rapid and accurate amplification, and optimizing the reaction composition, the success rate of target amplification in high-speed reactions has been improved while maintaining the fidelity of PrimeSTAR Max DNA Polymerase (Cat. #R045A/B). PrimeSTAR Max DNA Polymerase Ver.2 is a fully compatible and superior alternative to the existing one. Since this product is a 2X premixed reagent with an antibody-mediated hot start formulation to suppress DNA Polymerase activity and 3'  $\rightarrow$  5' exonuclease activity, it allows quick preparation of reactions at room temperature and is useful for high-throughput applications.

## <Features>

- Ultimate fidelity
- High-speed reaction (extension)
  - \* It can handle reactions involving large amounts of template DNA, which would ordinarily be difficult to amplify in high-speed reactions.
- High yield to the existing one
- Convenient premixed reagent, allowing easy reaction preps at room temperature

**Storage:** -20°C

**Note:** Gently mix well before use and centrifuge briefly. Repeated freeze-thaw of the Premix may reduce its activity.

## Electrophoresis, Cloning, and Sequencing of Amplified Products:

### (1) Electrophoresis

TAE Buffer is recommended for agarose gel electrophoresis of amplified products that are obtained using PrimeSTAR Max DNA Polymerase Ver.2.

**Note:** Use of TBE Buffer may result in DNA band patterns that are enlarged at the bottom of the gel.

### (2) Termini of amplified products

Most PCR products amplified with PrimeSTAR Max DNA Polymerase Ver.2 have blunt-end termini. Accordingly, they can be cloned directly into blunt-end vectors. If necessary, phosphorylate the amplified products before cloning. Use of Mighty Cloning Reagent Set (Blunt End) (Cat. #6027) is recommended for cloning into a blunt-end vector. In addition, In-Fusion cloning system is also applicable.

### (3) Restriction enzyme reaction

Prior to performing restriction enzyme digestion of amplified PCR products, remove all traces of PrimeSTAR Max DNA Polymerase Ver.2 from the reaction mixture by phenol/chloroform extraction or by using NucleoSpin Gel and PCR Clean-up (Cat. #740609.10/.50/.250)\*. Particularly for 3'-protruding restriction enzymes such as *Pst* I, the 3'-protruding termini produced by these enzymes may be deleted by 3'  $\rightarrow$  5' exonuclease activity of PrimeSTAR Max DNA Polymerase Ver.2, if residual polymerase remains present in the restriction digest reaction.

### (4) Direct sequencing

Perform phenol/chloroform extraction of PCR products prior to direct sequencing to ensure inactivation of 3'  $\rightarrow$  5' exonuclease activity. Alternatively, NucleoSpin Gel and PCR Clean-up\* may be used to purify DNA prior to sequencing.

\* Not available in all geographic locations. Check for availability in your area.

## General Composition of PCR Reaction Mixture (Total 50 $\mu$ l):

PrimeSTAR Max Ver.2 Premix (2X)	25 $\mu$ l (final conc. 1X)
Primer 1	10 - 15 pmol (final conc. 0.2 - 0.3 $\mu$ M)
Primer 2	10 - 15 pmol (final conc. 0.2 - 0.3 $\mu$ M)
Template	See "Template" chapter
Sterile purified water	up to 50 $\mu$ l

**Note:** The PCR reaction mixture can be prepared at room temperature. However, keep each of the reaction components on ice during the preparation process.

## PCR Conditions:

When performing rapid amplification protocols using PrimeSTAR Max DNA Polymerase Ver.2, 3-step reactions are recommended for best results.

For amplification of long or GC-rich targets, 2-step reactions are recommended when enhanced specificity is desired.

98°C *1	10 sec	} 30 - 35 cycle
55°C	5 or 15 sec*2	
68°C *3	5 - 10 sec/kb*4	
or		
98°C *1	10 sec	} 30 - 35 cycle
68°C *3	5 - 10 sec/kb*4	

### \*1 Denaturing conditions

During cycling, denaturation at 98°C for 5 to 10 sec is recommended. Denaturation at 94°C is also possible, but the time should be extended to 10 to 15 sec.

### \*2 Annealing time

For primer T<sub>m</sub> values of 55°C or greater, anneal for 5 sec.  
For primer T<sub>m</sub> values less than 55°C, anneal for 15 sec.

### [Important note]

Because the priming efficiency of PrimeSTAR Max DNA Polymerase Ver.2 is extremely high, longer annealing times may cause smearing of PCR products visible during electrophoresis analysis.

### \* T<sub>m</sub> value calculation:

$$T_m (^{\circ}\text{C}) = 2(\text{NA} + \text{NT}) + 4(\text{NC} + \text{NG}) - 5$$

where N represents the number of primer nucleotides having the specified identity (A, T, C, or G)

### \*3 Extension temperature

Set to 68°C for both 2-step and 3-step PCR.

### \*4 Extension time (as a guideline)

Amplified Product Sizes  $\leq$  6 kb : 5 sec/kb (10 sec/kb for cDNA)  
Amplified Product Sizes > 6 kb : 10 sec/kb (20 sec/kb for cDNA)  
In some cases, increasing the extension time can increase the amplification yield.

\* Select primer sequences using primer design software such as OLIGO Primer Analysis Software (Molecular Biology Insights, Inc.).

\* Do not use inosine-containing primers with PrimeSTAR Max DNA Polymerase Ver.2.

## Template:

Recommended quantities of template DNA (for 50  $\mu$ l reaction)

Human genomic DNA	5 - 500 ng
<i>E. coli</i> genomic DNA	100 pg - 200 ng
Plasmid DNA	10 pg - 10 ng
cDNA	Equivalent of 10 - 1,000 ng of total RNA

\* Do not use templates containing uracil, such as bisulfite-treated DNA.

PrimeSTAR is a registered trademark of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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