TaKaRa Ex Taq™ Hot Start Version

Code No. RR006A
Size : 250 units

Shipping at – 20°C
Store at – 20°C

Lot No.
Conc. : units/μl
Volume : μl
Expiration Date :

Description : *TaKaRa Ex Taq™* HS is designed to be suitable for Hot Start PCR. It is derived from *TaKaRa Ex Taq™* and neutralizing monoclonal antibody to *Taq* DNA polymerase. Non-specific amplification due to mis-priming and/or formation of primer dimer before thermal cycling can be prevented, since the antibody inhibits the polymerase activity by binding to the *Taq* DNA polymerase until the temperature elevates. This enzyme can be also used in general PCR conditions, since monoclonal antibody is denatured in the initial DNA-denaturation step.

Components : *TaKaRa Ex Taq™* HS (5 units/μl) 50 μl
10 × *Ex Taq* Buffer (20 mM Mg²⁺ plus) 1 ml
dNTP Mixture (2.5 mM each) 800 μl

Storage Buffer : 20 mM Tris-HCl (pH8.0)
100 mM KCl
0.1 mM EDTA
1 mM DTT
0.5% Tween 20
0.5% NP-40
50% Glycerol

Unit definition : One unit is the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition:

25 mM TAPS (pH 9.3 at 25°C)
50 mM KCl
2 mM MgCl₂
0.1 mM DTT
200 μM each dATP·dGTP·dCTP
100 μM [³H]-dTTP
0.25 mg/ml activated salmon sperm DNA

Purity :

Nicking, endonuclease and exonuclease activity were not detected after incubation of 0.6 μg of supercoiled pBR322 DNA, 0.6 μg of λ DNA or 0.6 μg of λ-*Hind* III digest with 10 units of this enzyme for 1 hour at 74°C.

Application : For DNA amplification by Hot Start PCR.

PCR products : As most PCR products amplified with *TaKaRa Ex Taq™* HS have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

PCR test : Good performance of DNA amplification by PCR was confirmed by using λ DNA (amplified fragment: 20 kb) and human genomic DNA (amplified fragment: β-globin gene 17.5 kb) as the templates. Inhibition of *Taq* DNA polymerase activity by the antibody was confirmed to be more than 90% after the reaction at 55°C for 10 min.

General reaction mixture for PCR (total 50 μl)

TaKaRa Ex Taq™ HS (5 units/μl) 0.25 μl
10 × *Ex Taq* Buffer 5 μl
dNTP Mixture (2.5 mM each) 4 μl
Template < 500 ng
Primer 1 10 - 50 pmol (final conc. 0.2 - 1.0 μM)
Primer 2 10 - 50 pmol (final conc. 0.2 - 1.0 μM)
Sterilized distilled water up to 50 μl

(Note) Reaction mixture can be set up at room temperature. Be sure to keep all reagents on ice.

PCR conditions :

This enzyme can be used in general PCR conditions, since the monoclonal antibody is denatured in the initial DNA-denaturation step. No need for a special step to denature the antibody to *Taq* polymerase.

Example : Amplification of 1 kb DNA fragment

98°C 10 sec. 98°C, 10 sec.] 30 cycles
55°C 30 sec. 30 cycles or 68°C, 1 min.]
72°C 1 min.

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. The recommendation is for 5 - 10 sec. at 98°C, or 20 - 30 sec. at 94°C.

NOTICE TO PURCHASER: LIMITED LICENSE

[P1] PCR Notice

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,079,352, 5,789,224, 5,618,711, 6,127,155 and claims outside the US corresponding to US Patent No. 4,889,818. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim (such as the patented 5' Nuclease Process claims in US Patents Nos. 5,210,015 and 5,487,972), no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

[L15] Hot Start PCR

Licensed under U.S. Patent No. 5,338,671 and 5,587,287, and corresponding patents in other countries.

[M57] LA Technology

This product is covered by the claims 6-16 of U.S. Patent No. 5,436,149 and its foreign counterpart patent claims.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com.

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