TaKaRa LA Taq[®] with GC Buffer

Code No. RR02AG Size: 125 U
 Conc.: 5 U/ μ l

Supplied Reagents:

2X GC Buffer I (5 mM Mg²⁺ plus) 1.25 ml
2X GC Buffer II (5 mM Mg²⁺ plus) 1.25 ml
dNTP Mixture (2.5 mM each) 400 μ l

Storage Buffer :

20 mM Tris-HCl, pH 8.0
100 mM KCl
0.1 mM EDTA
1 mM DTT
0.5% Tween 20
0.5% Nonidet P-40
50% Glycerol

Storage: -20°C

Unit definition :

One unit is the amount of enzyme that will incorporate 10 nmol of dNTPs into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition :

25 mM TAPS (pH 9.3 at 25°C)
50 mM KCl
2 mM MgCl₂
0.1 mM DTT
200 μ M ea. dATP·dGTP·dCTP
100 μ M [³H]-dTTP
0.25 mg/ml activated salmon sperm DNA

Purity :

Nicking activity was not detected after the incubation of 1 μ g of supercoiled pBR322 DNA with 25 units of this enzyme for 1 hour at 74°C.

Endonuclease and exonuclease activity were not detected after the incubation of 1 μ g of λ DNA or λ -Hind III digest with 25 units of this enzyme for 16 hours at 74°C.

PCR products :

As most PCR products amplified with *TaKaRa LA Taq* have one A at the 3'-termini, the obtained PCR products can be directly cloned into T-vectors. But when cloning long products (> 5 kb) into T-vectors, the cloning efficiency decreases. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

Applications :

For DNA amplification by PCR. Specifically, for amplification of targets with high GC-content or that contain repetitive sequences.

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

General reaction mixture for PCR (50 μ l reaction volume) :

| | |
|--------------------------------------|-----------------------------|
| <i>TaKaRa LA Taq</i> (5 U / μ l) | 0.5 μ l |
| 2X GC Buffer I or II * | 25 μ l |
| dNTP Mixture (2.5 mM each) | 8 μ l |
| Template | < 1 μ g |
| Primer 1 | final conc. 0.2 - 1 μ M |
| Primer 2 | final conc. 0.2 - 1 μ M |
| Sterile purified water | up to 50 μ l |

* This product is supplied with two reaction buffers. First, use 2X GC Buffer I. If DNA amplification is not obtained, amplification may be improved by using 2X GC Buffer II.

PCR conditions :

Amplification of a 1,255 bp region from human genomic DNA (total GC content: 65%)

| | | |
|------|--------|-------------|
| 94°C | 1 min | } 30 cycles |
| ↓ | | |
| 94°C | 30 sec | |
| 60°C | 30 sec | |
| 72°C | 2 min | } |
| ↓ | | |
| 72°C | 5 min | |

< Cool Start Method >

The "Cool Start Method" provides more accurate amplification and nonspecific amplification. This is a simple method that does not require specialized enzymes or additional reagents.

Cool Start Method Protocol

- 1) Keep all reagents on ice until use.
- 2) Prepare the reaction mixture on ice. *1,2
 - * 1 Order of reagent addition does not influence results.
 - * 2 Results will not be affected by leaving the mixture on ice for 30 min before thermal cycling.
- 3) Set a thermal cycler with the designated program. *3
 - * 3 PCR conditions do not need to be changed for Cool Start.
- 4) Set the tubes in a thermal cycler and start thermal cycling immediately.

TaKaRa LA Taq is a registered trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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