2XPCR Solution
Premix Ex Taq™ Hot Start Version

Code No. RR030A

Size: 500 µl X 5

(for 100 PCR reactions)

Shipping at -20°C

Stored at -20°C

Storage: Repeated freeze-thaw cycles may decrease the enzyme activity. Once it thawed, dispense into PCR tubes and store at -20°C. (ex. For 50 µl PCR reaction, dispense by 25 µl each tube.)

Lot No.

Expiry Date :

Description:

This product is an optimized mixture composed of enzyme (*TaKaRa Ex Taq™* HS), reaction buffer and dNTP mixture as 2-fold concentrations. *TaKaRa Ex Taq™* HS is designed to be suitable for Hot Start PCR. It is derived from *TaKaRa Ex Taq™* and neutralizing monoclonal antibody to *Taq* DNA polymerase. Non-specific amplification due to mispriming and/or formation of primer dimer before thermal cycling can be prevented, since the antibody inhibits the polymerase activity by binding to the *Taq* DNA polymerase until the temperature elevates. This enzyme can be used in general PCR conditions, since monoclonal antibody is denatured in the initial DNA-denaturation step.

Content:

TaKaRa Ex Taq™ HS *: 1.25 units/ 25 µl
dNTP Mixture : 2X conc.; ea. 0.4 mM
Ex Taq™ buffer : 2X conc.; including 4 mM Mg²⁺

***Specification of *TaKaRa Ex Taq™*HS (TaKaRa Cat.#RR006)**

Unit definition: One unit is the amount of the enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition:

25 mM TAPS (pH 9.3 at 25°C)
50 mM KCl
2 mM MgCl₂
1 mM 2-mercaptoethanol
200 µM each dATP, dGTP, dTTP
100 µM [α-³²P]-dCTP
0.25 mg/ml activated salmon sperm DNA

Purity: Nicking activity, endonuclease and exonuclease activity were not detected after the incubation of 0.6 µg of supercoiled pBR322 DNA, 0.6 µg of λDNA or 0.6 µg of λ-*Hind* III digest with 10 units of this enzyme for 1 hour at 74°C.

Test for antibody: Inhibition of *Ex Taq* activity by the antibody is confirmed to be more than 90% after the reaction at 55°C for 10 min.

Applications:

For DNA amplification by Polymerase Chain Reaction (PCR).

PCR product: As most PCR products amplified with *TaKaRa Ex Taq™* HS have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

PCR test : Good performance of DNA amplification by PCR was confirmed by using λDNA as the template (amplified fragment : 20 kbp). Good performance of DNA amplification of a single copy gene by PCR was also confirmed by using human genome DNA (amplified fragment : 2.9 kbp).

General reaction mixture for PCR (total 50 µl)

<i>Premix Ex Taq™</i> Hot Start Version*	25 µl
Template	<1 µg
Primer 1	0.2 ~ 1.0 µM (final conc.)
Primer 2	0.2 ~ 1.0 µM (final conc.)
Sterilized distilled water	up to 50 µl

* Please mix gently to be uniform and then use.

PCR conditions

This enzyme can be used in general PCR conditions, since the monoclonal antibody is denatured in the initial DNA-denaturation step. No need for a special step to denature the antibody to *Taq* polymerase.

(Example) Amplification of 1 kbp DNA fragment

98°C 10 sec.] 30 cycles	or	98°C 10 sec.] 30 cycles
55°C 30 sec.		68°C 1 min.		
72°C 1 min.				

Note:

Denaturation condition varies depending on an used thermal cycler and tube. It is recommended for 20 - 30 sec. at 94°C, or 5 - 10 sec. at 98°C.

Note

For research use only. Not for use in diagnostic or therapeutic procedures.

U.S. Patent 5,436,149 for LA Technology is owned by TAKARA BIO INC.

NOTICE TO PURCHASER: LIMITED LICENSE

A license under the foreign counterparts of U.S. Patents Nos. 4,683,202, 4,683,195 and 4,965,188 owned by F.Hoffmann-La Roche Ltd. and U.S. Patent No. 5,075,216 and its foreign counterpart, owned by Roche Molecular Systems, Inc. and F.Hoffmann-La Roche Ltd. for use in research and development, has an up-front fee component and a running-royalty component. The purchase price of this product includes limited, nontransferable rights under the running-royalty component to use only this amount of the product to practice the polymerase chain reaction (PCR) and related processes described in said patents where such processes are covered by patents solely for the research and development activities of the purchaser when this product is used in conjunction with a thermal cycler whose use is covered by the up-front fee component. Rights to the up-front fee component must be obtained by the end user in order to have a complete license to use this product in the PCR process where the process is covered by patents. These rights under the up-front fee component may be purchased from Applied Biosystems or obtained by purchasing an authorized thermal cycler. No right to perform or offer commercial services of any kind using PCR, where the process is covered by patents, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted by implication or estoppel. Further information on purchasing licenses to practice the PCR process where the process is covered by patents may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or the Licensing Department, Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

2 × PCR Solution
Premix Ex Taq™ Hot Start Version

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