# TaKaRa LA Taq® Hot Start Version

**Code No. RR042A**      **Size:**            **125 U**  
                                 **Conc.:**            **5 U/μl**

## Supplied Reagents:

**10X LA PCR Buffer II (Mg<sup>2+</sup> plus) (25 mM)**    **1 ml**  
**dNTP Mixture (2.5 mM each)**                    **400 μl**

## Description:

*TaKaRa LA Taq* HS is designed for hot start PCR; it includes a neutralizing monoclonal antibody that recognizes *LA Taq* DNA polymerase. This antibody inhibits polymerase activity by binding to *LA Taq*, thereby preventing nonspecific amplification due to mispriming and/or formation of primer dimers during reaction set-up before thermal cycling. Antibody-mediated repression is released during the initial DNA denaturation step of PCR. This enzyme can be used with standard PCR conditions.

## Storage Buffer:

20 mM	Tris-HCl, pH 8.0
100 mM	KCl
0.1 mM	EDTA
1 mM	DTT
0.5%	Tween 20
0.5%	NP-40
50%	Glycerol

**Storage:**            -20°C

## dNTP Mixture:

dNTP Mixture is ready for use in PCR without dilution.

Concentration	: 2.5 mM of each dNTP
Form	: Dissolved in water (sodium salts), pH 7 - 9
Purity	: ≥ 98% for each dNTP

## Unit definition:

One unit is the amount of enzyme that will incorporate 10 nmol of dNTPs into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

## Reaction mixture for unit definition:

25 mM	TAPS (pH 9.3 at 25°C)
50 mM	KCl
2 mM	MgCl <sub>2</sub>
0.1 mM	DTT
200 μM	each dATP·dGTP·dCTP
100 μM	[ <sup>3</sup> H]-dTTP
0.25 mg/ml	activated salmon sperm DNA

## Purity:

Nicking, endonuclease, and exonuclease activity were not detected after the incubation of 0.6 μg of supercoiled pBR322 DNA, 0.6 μg of λ DNA, or 0.6 μg of λ-*Hind* III digest with 10 units of this enzyme for 1 hour at 74°C.

## Application:

For DNA amplification by hot start PCR. Specifically, for amplifying fragments > 15 kb.

## PCR products:

As most PCR products amplified with *TaKaRa LA Taq* HS have one A at the 3'-termini, the obtained PCR products can be directly cloned into T-vectors. When cloning long products (> 5 kb) into T-vectors, the cloning efficiency may be very low. It is also possible to clone the product in blunt-end vectors after blunting and phosphorylation of the ends.

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## General reaction mixture for PCR (50 μl reaction volume):

<i>TaKaRa LA Taq</i> HS (5 U/μl)	0.5 μl
10X LA PCR Buffer II (Mg <sup>2+</sup> plus)	5 μl
dNTP Mixture (2.5 mM each)	8 μl
Template	< 1 μg
Primer 1	0.2 - 1.0 μM (final conc.)
Primer 2	0.2 - 1.0 μM (final conc.)
Sterile purified water	up to 50 μl

## PCR conditions:

This enzyme can be used with standard PCR conditions, since the monoclonal antibody is denatured in the initial DNA-denaturation step. There is no need for an additional step to denature the anti-*LA Taq* antibody.

Example: Amplification of a 17.5 kb DNA fragment

94°C	1 min	} 30 cycles
↓		
98°C	10 sec	
68°C	15 min	}
↓		
72°C	10 min	

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. Denaturation for 5 - 10 sec at 98°C, or 20 - 30 sec at 94°C is recommended.

*TaKaRa LA Taq* is a registered trademark of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takara-bio.com](http://www.takara-bio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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