

For Research Use

Takara

One Step TB Green[®] PrimeScript[™] RT-PCR Kit (Perfect Real Time)

Product Manual

We have begun the process of changing the names for Takara Bio's intercalator-based real-time PCR (qPCR) products to the "TB Green series". These products can be used the same way as before, as only the names are changing. Catalog number and product performance are unaffected by this transition.

Table of Contents

| | | |
|-------|-----------------------------------|----|
| I. | Description..... | 3 |
| II. | Principle | 3 |
| III. | Components | 5 |
| IV. | Storage | 6 |
| V. | Features..... | 6 |
| VI. | Precautions..... | 6 |
| VII. | Protocol..... | 7 |
| VIII. | Optimization | 12 |
| IX. | Experimental Example | 13 |
| X. | Appendix..... | 14 |
| XI. | Guidelines for Primer Design..... | 15 |
| XII. | Related Products | 16 |

I. Description

The One Step TB Green PrimeScript RT-PCR Kit (Perfect Real Time) is designed for intercalator-based one step, real-time RT-PCR using TB Green for detection. All reactions can be performed in a single tube, therefore operation is simple, and the risk of contamination is minimized.

Amplified products are monitored in real time, so there is no need for electrophoresis after PCR. This kit is suitable for detection of tiny amounts of RNA such as RNA virus.

This kit uses PrimeScript RTase, which has excellent extension capability and can efficiently synthesize cDNA in short time, and *Takara Ex Taq*® HS, a high efficiency hot start PCR enzyme. Both enzymes are optimized for one step RT-PCR, allowing efficient generation of RT-PCR products.

Applicable real time PCR instruments:

- Thermal Cycler Dice™ Real Time System III (Cat. #TP950/TP970/TP980/TP990)*
- Thermal Cycler Dice Real Time System II (Cat. #TP900/TP960)*
- Thermal Cycler Dice Real Time System *Lite* (Cat. #TP700/TP760)*
- Applied Biosystems 7300/7500/7500 Fast Real-Time PCR System, StepOnePlus Real-Time PCR System (Thermo Fisher Scientific)
- Smart Cycler System/Smart Cycler II System (Cepheid)
- LightCycler (Roche Diagnostics), and other

* Not available in all geographic locations. Check for availability in your area.

II. Principle


The One Step TB Green PrimeScript RT-PCR Kit (Perfect Real Time) enables cDNA synthesis from RNA using reverse transcriptase, PrimeScript RTase, and PCR amplification using *Takara Ex Taq* HS in one tube. PCR amplification products are monitored in real time with TB Green as an intercalator

1. PCR

PCR is a technique that amplifies a targeted DNA region from a small amount of DNA. One cycle of PCR includes heat denaturation of DNA, primer annealing, and primer extension by DNA polymerase. By repeating this process, PCR allows amplification of a targeted gene segment to one million times in a short time period. The use of a hot start PCR enzyme, *Takara Ex Taq* HS, for amplification prevents miss-priming during preparation of reagents, and non-specific amplification caused by primer dimer.

2. RT-PCR

Although RNA cannot be used as a direct template of PCR, PCR can be applied for RNA analysis when cDNA is first synthesized from RNA with reverse transcriptase (RT-PCR). The principle of one step RT-PCR is shown below; a specific primer (reverse) is used for reverse transcription. Then using the synthesized cDNA as a template, PCR amplification is performed with specific primers (forward, reverse). (Random Primer and Oligo dT Primer cannot be used for reverse transcription.)




Principle of One Step RT-PCR

3. Fluorescence detection


[Intercalator method]

This method detects fluorescence produced during amplification by adding an Intercalator (e.g., TB Green) that fluoresces when bound to double strand of DNA. By detecting this fluorescence, the amount of amplified DNA can be quantified and the melting point can be measured.


1) Heat denaturation



2) Primer annealing



3) Extension



III. Components (100 reactions, 50 µl volume)

| | |
|---|-------------|
| 1. 2X One Step TB Green RT-PCR Buffer III*1 | 840 µl x 3 |
| 2. <i>Takara Ex Taq</i> HS (5 U/µl) | 100 µl |
| 3. PrimeScript RT enzyme Mix II*2 | 100 µl |
| 4. RNase Free dH ₂ O | 1.25 ml x 2 |
| 5. ROX Reference Dye (50X conc.)*3 | 100 µl |
| 6. ROX Reference Dye II (50X conc.)*3 | 100 µl |

*1 Includes dNTP Mixture, Mg²⁺, and TB Green. Store protected from light.

*2 Includes PrimeScript RTase and RNase Inhibitor.

*3 ROX Reference Dye/Dye II is used for normalization of the fluorescent signal by background subtraction.

◆ Use the ROX Reference Dye

- Applied Biosystems 7300 Real-Time PCR System (Thermo Fisher Scientific)
- StepOnePlus Real-Time PCR System (Thermo Fisher Scientific)

◆ Use the ROX Reference Dye II

- Applied Biosystems 7500/7500 Fast Real-Time PCR System (Thermo Fisher Scientific)

◆ Not required

- Thermal Cycler Dice Real Time System series (Cat. #TP950/TP900/TP700 etc.)
- Smart Cycler System/Smart Cycler II System (Cepheid)
- LightCycler (Roche Diagnostics)

Materials required but not provided

1. Real-time PCR instrument
2. Reaction tube or plates suitable for the real-time PCR instrument
3. PCR primers*
4. Micropipettes and pipette tips (autoclaved)

* Refer to XI. Guidelines for Primer Design.

IV. Storage

-20°C

Note: Store 2X One Step TB Green RT-PCR Buffer III protected from light.

V. Features

- (1) One Step RT-PCR is intended for accurate and rapid analysis of RNA viruses or small amounts of RNA.
- (2) For PCR, *TaKaRa Ex Taq* HS, a high efficiency hot start PCR enzyme, is used. The buffer system is optimized for real-time PCR, thus high quality detection is possible. Moreover, One Step TB Green RT-PCR Buffer III is a 2X concentration of premix with TB Green, so preparation of reaction mixtures is simple.

VI. Precautions

Please read it carefully before use.

- (1) When mixing reagents for PCR, mix enough for 10 reactions (master mix). Using a master mix allows accurate reagent dispensing, minimizes reagent pipetting errors, and avoids repeated dispensing of each reagent. This helps to minimize variation between experiments.
- (2) PrimeScript RT enzyme Mix II and *TaKaRa Ex Taq* HS should be mixed gently. Avoid generating bubbles. Gently spin down the solution prior to pipetting. Pipet the enzymes slowly as the enzyme contains 50% glycerol and is very viscous. Keep the enzyme at -20°C until just before use and return into the freezer promptly after use.
- (3) If precipitate appears during thawing of 2X One Step TB Green RT-PCR Buffer III, dissolve completely by vortex before using.
- (4) Use new disposable pipette tips to avoid contamination between samples for transferring reagent.
- (5) Use the gene-specific reverse primer for reverse transcription. Random Primer or Oligo-dT Primer should not be used.

VII. Protocol

[For the Thermal Cycler Dice Real Time System III, II, and Lite]

1. Prepare the following reagents on ice.

< Per reaction >

| Reagent | Volume | Final conc. |
|--|-----------------------------|---------------|
| 2X One Step TB Green RT-PCR Buffer III | 12.5 μ l | 1X |
| Takara Ex Taq HS (5 U/ μ l) | 0.5 μ l | |
| PrimeScript RT enzyme Mix II | 0.5 μ l | |
| PCR Forward Primer (10 μ M) | 0.5 μ l | 0.2 μ M*1 |
| PCR Reverse Primer (10 μ M) | 0.5 μ l | 0.2 μ M*1 |
| total RNA | 2 μ l | *2 |
| RNase Free dH ₂ O | 8.5 μ l | |
| Total | 25 μl | |


*1 The final concentration of primers can be 0.2 μ M in most reactions. However, if there is insufficient reactivity, use a primer concentration between 0.1 and 1.0 μ M.

*2 It is recommended to use 10 pg - 100 ng total RNA as a template.

2. Gently spin down the reaction tubes or plate, then start the reaction after setting them in the Thermal Cycler Dice Real Time System.

Use the standard protocol described as follows. Try this protocol first and optimize PCR condition as necessary. Perform 3-step PCR when shuttle PCR is difficult, for example for primer with low T_m value.

To optimize PCR conditions further, refer to section VIII. Optimization.



Pattern 1: Reverse Transcription

Hold
42°C 5 min
95°C 10 sec

Pattern 2: PCR

Cycles: 40
95°C 5 sec
60°C 30 sec

Pattern 3: Dissociation

Note:

This product contains *Takara Ex Taq HS*, which includes a *Taq* antibody. There is no need to heat at 95°C for 5 - 15 min for initial denaturation (as is required for chemically modified *Taq* polymerases). If longer heat treatment is used with this kit, the enzyme activity decreases and amplification efficiency and the accuracy for quantification may also be affected. Perform heat inactivation of RTase prior to PCR at 95°C for 10 sec.

3. After the reaction is complete, verify the amplification curve and melting curve. Plot a standard curve if absolute quantification is desired. Refer to the instrument manual of Thermal Cycler Dice Real Time System and the following application examples for the analysis method with Thermal Cycler Dice Real Time System.

[StepOnePlus Real-Time PCR System]

* Follow manuals of each instrument to operate them.

1. Prepare the following reagents on ice.

< Per reaction >

| Reagent | Volume | Volume | Final Conc. |
|--|-------------------------------|-------------------------------|---------------|
| 2X One Step TB Green RT-PCR Buffer III | 10 μ l | 25 μ l | 1X |
| TaKaRa Ex Taq HS (5 U/ μ l) | 0.4 μ l | 1 μ l | |
| PrimeScript RT enzyme Mix II | 0.4 μ l | 1 μ l | |
| PCR Forward Primer (10 μ M) | 0.4 μ l | 1 μ l | 0.2 μ M*1 |
| PCR Reverse Primer (10 μ M) | 0.4 μ l | 1 μ l | 0.2 μ M*1 |
| ROX Reference Dye or Dye II (50 X)*3 | 0.4 μ l | 1 μ l | |
| total RNA | 2 μ l | 4 μ l | *2 |
| RNase Free dH ₂ O | 6 μ l | 16 μ l | |
| Total | 20 μl*4 | 50 μl*4 | |

*1 The final concentration of primers can be 0.2 μ M in most reactions. However, if there is insufficient reactivity, use a primer concentration between 0.1 and 1.0 μ M.

*2 It is recommended to use 20 pg - 200 ng total RNA as a template in a 50 μ l reaction.

*3 The ROX Reference Dye/Dye II is supplied for performing normalization of fluorescent signal intensities between wells when using the real-time PCR instruments below. For Applied Biosystems 7300 Real-Time PCR System and StepOnePlus Real-Time PCR System, use ROX Reference Dye (50X). For Applied Biosystems 7500/7500 Fast Real-Time PCR System, use ROX Reference Dye II (50X).


*4 Prepare in accordance with the recommended volume for each instrument.

2. Start reaction

Use the standard protocol described below. Try this protocol first and optimize PCR conditions as necessary. Perform 3-step PCR when shuttle PCR is difficult, for example for primer with low T_m value.

To optimize PCR conditions further, refer to section VIII. Optimization.

<Applied Biosystems 7300/7500 Real-Time PCR System and StepOnePlus Real-Time PCR System>



Stage 1, 2: Reverse transcription

Reps: 1
42°C 5 min
95°C 10 sec

Stage 3: PCR

Reps: 40
95°C 5 sec
60°C 30 - 34 sec*

Stage 4: Dissociation Protocol

* Set up at 30 sec for StepOnePlus, 31 sec for 7300, and 34 sec for 7500.

< Applied Biosystems 7500 Fast Real-Time PCR System >

Holding Stage

42°C 5 min

95°C 10 sec

Cycling Stage

Number of Cycles: 40

95°C 3 sec

60°C 30 sec

Melt Curve Stage

Note:

This product contains *TAKARA Ex Taq* HS, which includes a *Taq* antibody. There is no need to heat at 95°C for 5 - 15 min for initial denaturation (as is required for chemically modified *Taq* polymerases). If longer heat treatment is used with this kit, the enzyme activity decreases and amplification efficiency and the accuracy for quantification may also be affected. Perform heat inactivation of RTase prior to PCR at 95°C for 10 sec.

3. After the reaction is complete, verify the amplification curve and dissociation curve. Plot the standard curve if quantification is done. Refer to the operation manual of the real-time PCR instrument for additional analysis methods.

[For the LightCycler]

* Follow the manual of LightCycler (Roche Diagnostics).

1. Prepare the following reagents on ice.

< Per reaction >

| Reagent | Volume | Final Conc. |
|--|-----------------------------|---------------|
| 2X One Step TB Green RT-PCR Buffer III | 10 μ l | 1X |
| Takara Ex Taq HS (5 U/ μ l) | 0.4 μ l | |
| PrimeScript RT enzyme Mix II | 0.4 μ l | |
| PCR Forward Primer (10 μ M) | 0.4 μ l | 0.2 μ M*1 |
| PCR Reverse Primer (10 μ M) | 0.4 μ l | 0.2 μ M*1 |
| total RNA | 2 μ l | *2 |
| RNase Free dH ₂ O | 6.4 μ l | |
| Total | 20 μl | |


*1 The final concentration of primers can be 0.2 μ M in most reactions. However, if there is insufficient reactivity, use a primer concentration between 0.1 and 1.0 μ M.

*2 It is recommended to use 10 pg - 100 ng total RNA as a template.

2. Gently spin down the PCR capillaries, then start the reaction after setting them onto LightCycler.

Use standard protocol described below. Try this protocol first and optimize PCR conditions as necessary. Perform 3-step PCR when shuttle PCR is difficult, for example for primers with low T_m values.

To optimize PCR conditions further, refer to section VIII. Optimization.



Stage 1: Reverse transcription
 42°C 5 min 20°C/sec
 95°C 10 sec 20°C/sec
 1 cycle
 Stage 2: PCR
 95°C 5 sec 20°C/sec
 60°C 20 sec 20°C/sec
 40 cycles
 Stage 3: Melting curve analysis
 95°C 0 sec 20°C/sec
 65°C 15 sec 20°C/sec
 95°C 0 sec 0.1°C/sec

Note:

This product contains *Takara Ex Taq HS*, which includes a *Taq* antibody. There is no need to heat at 95°C for 5 - 15 min for initial denaturation (as is required for chemically modified *Taq* polymerases). If longer heat treatment is used with this kit, the enzyme activity decreases and amplification efficiency and the accuracy for quantification may also be affected. Perform heat inactivation of RTase prior to PCR at 95°C for 10 sec.

3. After the reaction is complete, verify amplification curve and melting curve. Plot the standard curve if absolute quantification is desired. Refer to the operation manual of the LightCycler for additional analysis methods.

[For the Smart Cycler II System]

1. Prepare the following reagents on ice.

< Per reaction >


| Reagent | Volume | Final Conc. |
|--|-----------------------------|---------------|
| 2X One Step TB Green RT-PCR Buffer III | 12.5 μ l | 1X |
| <i>Takara Ex Taq HS</i> (5U/ μ l) | 0.5 μ l | |
| PrimeScript RT enzyme Mix II | 0.5 μ l | |
| PCR Forward Primer (10 μ M) | 0.5 μ l | 0.2 μ M*1 |
| PCR Reverse Primer (10 μ M) | 0.5 μ l | 0.2 μ M*1 |
| total RNA | 2 μ l | *2 |
| RNase Free dH ₂ O | 8.5 μ l | |
| Total | 25 μl | |

- *1 The final concentration of primers can be 0.2 μ M in most reactions. However, if there is insufficient reactivity, use a primer concentration between 0.1 and 1.0 μ M.
- *2 It is recommended to use 10 pg - 100 ng total RNA as a template.

2. Gently spin down the reaction tubes with a Smart Cycler specific centrifuge, then start the reaction after setting them onto Smart Cycler.
Use standard protocol described below. First, try this protocol and optimize PCR condition as necessary. Perform 3-step PCR when shuttle PCR is difficult to process, for example for primer with low T_m value.

To optimize PCR conditions further, refer VIII. Optimization.

[For the Applied Biosystems 7300/7500/7500 Fast Real-Time PCR System and



- Stage 1: Reverse Transcription Hold
 - 42°C 5 min
 - 95°C 10 sec
- Stage 2: PCR
 - Repeats: 40 times
 - 95°C 5 sec
 - 60°C 20 sec
- Stage 3: Melt Curve

Note:

This product contains *Takara Ex Taq HS*, which includes a *Taq* antibody. There is no need to heat at 95°C for 5 - 15 min for initial denaturation (as is required for chemically modified *Taq* polymerases). If longer heat treatment is used with this kit, the enzyme activity decreases and amplification efficiency and the accuracy for quantification may also be affected.
Perform heat inactivation of RTase prior to PCR at 95°C for 10 sec.

3. After the reaction is complete, verify the amplification curve and melting curve. Plot the standard curve if absolute quantification is desired.
Refer to the operation manual of Smart Cycler for additional analysis methods.

VIII. Optimization


If the recommended conditions (shuttle PCR standard protocol) provide unsatisfactory reactivity, follow the procedures below to evaluate what primer concentration and PCR conditions to use.

Select PCR conditions based on comprehensive analysis, taking into consideration both reaction specificity and amplification efficiency. A PCR system balanced between these two aspects allows accurate assay over a wide range of concentrations.

- System with a high reaction specificity
 - With no template control, non-specific amplification such as primer-dimer formation do not take place.
 - Non-specific amplification products, those other than the target product, are not generated.
- System with a high amplification efficiency
 - Amplification product is detected at early cycles (small Ct value).
 - PCR amplification efficiency is high (near the theoretical value of 100%).

1. Evaluation of primer concentration


The relationship between primer concentration and reaction specificity and amplification efficiency is illustrated below. Reducing the primer concentration raises reaction specificity. Increasing the primer concentration, in contrast raises amplification efficiency.



2. Evaluation of PCR condition


- To raise reaction specificity

Raising the annealing temperature may improve reaction specificity. Perform optimization while checking amplification efficiency.



- To raise amplification efficiency

Prolonging the elongation time or switching to a 3-step PCR may improve amplification efficiency. Perform optimization using the steps below.




IX. Experimental Example

Detection of Rat Rplp2 (ribosomal protein, large P2)
(Thermal Cycler Dice Real Time System is used here.)


1. Procedure

Using total RNA 6.4 pg - 100 ng that was prepared from rat liver or sterile purified water (negative control) as a template, real-time one step RT- PCR was performed.

2. Result




Amplification Curve




Melting Curve

After reaction, obtain Ct value from 2nd derivative of amplification curve and plot the standard curve.



2nd derivative



Standard Curve

Target DNA was detected. Using 6.4 pg - 100 ng of total RNA, the melting curve shows that the same amplification products were obtained even different amounts of template. The linearity of standard curve was obtained over the wide range of template concentrations.

X. Appendix

Preparation of RNA sample

This kit is designed to perform the reverse transcription of RNA to cDNA and subsequent PCR amplification. It is important to use high purity of RNA sample for better yields of cDNA synthesis. So, it is essential to inhibit cellular RNase activity and also to prevent contamination with RNases derived from equipment and solutions used. Extra precaution should be taken during sample preparation, including use of clean disposable gloves, dedication of a work surface exclusively for RNA preparation, and avoiding unnecessary talking during assembly, to prevent the RNase contamination from operators' sweat or saliva.

[Equipment]

Disposable plastic equipment should be used. Glass tools should be treated with the following protocol prior to use.

- (1) Hot-air sterilization (180°C, 60 min)
- (2) Treatment with 0.1% DEPC at 37°C, for 12 hours followed by autoclaving at 120°C for 30 min to remove DEPC.

* It is recommended that all equipment are used exclusively for RNA preparation.

[Reagent]

All reagents used in this experiment must be prepared using tools treated as described above (hot-air sterilization(180°C, 60 min) or DEPC treatment), and purified water must be treated with 0.1% DEPC and autoclaved. All reagents and purified water should be used exclusively for RNA experiments.

[Preparation of RNA sample]

Since RT-PCR usually requires only small amounts of RNA, common purification methods are usually sufficient. However, we recommend that the guanidine thiocyanate (GTC) method be used if possible. In general, RNA should be of the highest purity possible.

When preparing high purity total RNA from cell cultures or tissue samples, NucleoSpin RNA (Cat. #740955.10/.50/.250) or the AGPC method simplified reagent RNAiso Plus (Cat. #9108/9109) can be used. For blood samples, NucleoSpin RNA Blood (Cat. #740200.10/.50) or RNAiso Blood (Cat. #9112/9113) can be used.

XI. Guidelines for Primer Design

It is essential to design primers that allow good reactivity for successful real-time PCR reactions. Follow the guidelines below to design primers which offer high amplification efficiency and minimize non-specific reaction.

Amplification product

| | |
|----------------|--|
| Amplified size | The optimal size is 80 - 150 bp (amplification up to 300 bp is possible) |
|----------------|--|

Primer

| | |
|--------------------|--|
| Length | 17 - 25 mer |
| GC content | 40 - 60% (45 - 55% is recommended.) |
| Tm | Make sure that the Tm values for the forward primer and the reverse primer do not differ greatly. Use primer design software to determine Tm values. OLIGO*1: 63 - 68°C Primer 3: 60 - 65°C |
| Sequence | Make sure that overall there are no base sequence biases. Avoid having any GC-rich or AT-rich regions in the sequence (particularly at the 3' end). Avoid having consecutive T/C pairings (polypyrimidine). Avoid having consecutive A/G pairings (polypurine). |
| Sequence of 3' end | Avoid having any GC-rich or AT-rich sequence at the 3' end. It is preferable to have a G or C as the 3' end-base. Avoid primers with T as the 3' end-base. |
| Complementarity | Avoid having any complementary sequences of 3 bases or more within a primer and between primers. Avoid having any complementary sequences of 2 bases or more at the primer's 3' ends. |
| Specificity | Verify primer specificity by BLAST search*2. |

*1 OLIGO Primer Analysis Software (Molecular Biology Insights, Inc.)

*2 <http://www.ncbi.nlm.nih.gov/BLAST/>

XII. Related Products

One Step TB Green® PrimeScript™ PLUS RT-PCR Kit (Perfect Real Time) (Cat. #RR096A/B)*
One Step TB Green® PrimeScript™ RT-PCR Kit II (Perfect Real Time) (Cat. #RR086A/B)
One Step PrimeScript™ RT-PCR Kit (Perfect Real Time) (Cat. #RR064A/B)*
TB Green® *Premix Ex Taq*™ (Tli RNaseH Plus) (Cat. #RR420/A/B)
TB Green® *Premix Ex Taq*™ II (Tli RNaseH Plus) (Cat. #RR820/A/B)
TB Green® Premix DimerEraser™ (Perfect Real Time) (Cat. #RR091A/B)*
Probe qPCR Mix (Cat. #RR391A/B)*
Thermal Cycler Dice™ Real Time System III (Cat. #TP950/TP970/TP980/TP990)*
Thermal Cycler Dice™ Real Time System II (Cat. #TP900/TP960)*
Thermal Cycler Dice™ Real Time System *Lite* (Cat. #TP700/TP760)*
RNAiso Plus (Cat. #9108/9109)
RNAiso Blood (Cat. #9112/9113)

* Not available in all geographic locations. Check for availability in your area.

TB Green and *Takara Ex Taq* are registered trademarks of Takara Bio Inc.
PrimeScript, Thermal Cycler Dice, DimerEraser, and *Premix Ex Taq* are trademarks of Takara Bio Inc.

NOTE: This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.
