Tris Buffer Powder, pH 8.0

Code No. T9152 **Size: 10 pouches**
(for 10 L)

* 2 years from date of receipt under proper storage conditions.

Description :

Tris buffer is used in many applications in biochemistry and molecular biology.
This product is a powder for preparing 1 M Tris Buffer (pH 8.0). The buffer can easily be prepared by dissolving it with H₂O.
One pouch is for preparing 1,000 ml of 1 M Tris Buffer (pH 8.0) .

Storage :

Room temperature at desiccated condition.

Specifications :

Format : Powder
Concentration : 1 M Tris-HCl
Volume : For 1,000 ml
pH : 8.0 ± 0.1 at 25°C

Preparation for use :

Transfer one pouch of Tris Buffer Powder, pH 8.0 in a flask or beaker.
Add 300 ml of distilled water or deionized water and stir the solution for a few minutes. Adjust the volume to 1,000 ml, and stir until complete dissolution.

Applications :

For preparation of various buffer (e.g. TE Buffer)

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.
Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.
If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.
Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.
All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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