

## Table of contents

1.	Description	2
2.	Principle	2
3.	Feature	2
4.	Kit component	2
5.	Storage	3
6.	Reference	3
7.	Available specimen form	3
8.	Preparation of labeling reaction mixture	3
9.	Reagents and instruments required other than this kit	3
10.	Protocol	
	A. Cultured cell	4
	B. Frozen tissue section	5
	C. Paraffin-embedded tissue section	5
	D. Deparaffinization	6
11.	Note	6
12.	Flow chart of procedure	7
13.	Control Slide Experiment	8
14.	Q&A	8

### *Also available from TaKaRa*

---

ApopLadder Ex™ .....Cat.#MK400



## In situ Apoptosis Detection Kit

### 5. Storage

Shipped at -20°C .

Store the component separately

1, 2, 5	.....at -20°C
3	.....at 4°C *
4	.....at room temperature**

\*Store at 4°C once it thawed.

\*\*Store at room temperature after delivered, though the kit is shipped at -20°C .

### 6. Reference:

1. Dawn, R.Z. Tornusciolo, Robert E.Schmidt, and Kevin A.Roth (1995) *BioTechniques*, **19**, 800-805.
2. Bernard Piqueras, Brigitte Autran, Patrice Debre, and Guy Gorochov (1996) *BioTechniques*, **20**, 634-640.
3. Tushar Patel, Amindra Arora, and Gregory J.Gores (1995) *Analytical Biochemistry*, **229**, 229-235.
4. Jan H.Wijsman, Richard R. Jonker, Rob Keijer, Cornelis J.H. Van De Velde, Cees J. Cornelisse, and Jan Hein Van Dierendonck (1993) *The Journal of Histochemistry and Cytochemistry*, **41**, No.1, 7-12.
5. R. Gold, M. Schmied, G. Rothe, H.Zischler, H. Breitschopf, H.Wekerle, and H. Lassmann (1993) *The Journal of Histochemistry and Cytochemistry*, **41**, No.7, 1023-1030.
6. Yael Garvrieli, Yoav Sherman, and Shmuel A. Ben-Sasson (1992) *The Journal of Cell Biology*, **119**, No.3, 493-501.
7. Jorn Strater, Andreas R, Gunthert, Silke Bruderlein, Peter Moller (1995) *Histochemistry* **103**, 157-160.
8. Ruth J. Muschel, Eric J. Bernhard, Luis Garza, W. Gillies McKenna, and Cameron J. Koch (1995) *Cancer Research*, **55**, 995-998.
9. Katsuaki Kato (1996) *Progress in Medicine*, **16**, No.3, 867-870.
10. Xun Li, Jianping Gong, Eric Feldman, Karen Seiter, Frand Traganos and Zbigniew Darzynkiewicz (1994) *Leukemia and Lymphoma*, **13**, Suppl. 1, 65-70.
11. Alexander Dolzhanskiy, Ross S. Basch (1995) *Journal of Immunological Methods*, **180**, 131-140.
12. Gold R., Schmied M. *et al.* (1994) *Laboratory Investigation* **71** No.2, 219-225.

### 7. Available specimen form:

- Cell: Adherent cells (Cultured on a chamber slide)  
Cells suspension (Cytospin preparation on a slide glass, or smeared and collected in a microtube)
- Tissue section: Frozen section, paraffin-embedded section

### 8. Preparation of labeling reaction mixture

- 1) Per one sample, prepare the reaction mixture by adding 5  $\mu$  l of TdT Enzyme to 45  $\mu$  l of Labeling Safe Buffer. One tube of Labeling Safe Buffer also contains buffer in the amount for 2 negative control reactions (50  $\mu$  l x 2).
- 2) Mix the prepared mixture gently but well to be mixed uniformly. Reaction mixture should be prepared each time use, and should be stored on ice until use. Do not store the prepared mixture. If it is left for a long period, the enzyme in the mixture might be inactivated.  
**Note:** Anti-FITC HRP Conjugate needs no preparation prior to use.

### 9. Reagents and instruments required other than this kit

Please refer to specific protocols to determine applicability.

- Distilled water
- Washing buffer (PBS or TBS)  
ex. Phosphate Buffered Salts (Cat.#T900)  
TBS (Tris-Buffered Saline) powder (Cat.#T903)
- Coloring substrate (DAB)  
(ex. SIGMA Chemical Co. (Code.#D-5905), DAKO(Code.#S3000))
- H<sub>2</sub>O<sub>2</sub>
- Counterstain solution (methyl green)
- Micropipette and microtube (autoclaved)
- Humidified chamber
- Incubator (37°C )
- Glass coverslip, Guard  
ex. Takara Slide Seal for in situ PCR (Cat.#9066, #9067, #9068)
- Cover glass
- Slide glass (precoated with silan)
- Microscope (Both fluorescence and light one are applicable.)

[For detection using paraffin-embedded section]

- Glass or plastic coplin jar
- Xylene
- Ethanol (100%, 90%, 80%)
- Glass coverslip  
ex. Takara Slide Seal for in situ PCR (Cat.#9066, #9067, #9068)
- Proteinase K (Cat.#9033) (20  $\mu$ g/ml)
- 3% H<sub>2</sub>O<sub>2</sub> (For endogenous peroxidase inactivation)
- Mounting medium

[For detection using frozen section]

- Slideglass pretreated to prevent exfoliation (precoated with silan)
- Fixation solution (ex. 10% neutral-buffered formalin, acetone, 4% paraformaldehyde, etc.)
- Methanol containing 0.3% H<sub>2</sub>O<sub>2</sub> (for Endogenous peroxidase inactivation)

[For detection using cell]

- Slideglass pretreated to prevent exfoliation (precoated with silan)
- Fixation solution (ex. 10% neutral-buffered formalin, 4% paraformaldehyde, etc.)
- Methanol containing 0.3% H<sub>2</sub>O<sub>2</sub> (for Endogenous peroxidase inactivation)
- Microcentrifuge (cytospin)

## [A. Cultured cell]

## 10.PROTOCOL

1. Wash the collected cells with PBS, and dry in air on a silanized slideglass. Fix the cells with 4% paraformaldehyde / PBS solution (pH7.4) by leaving at room temperature for 15-30 min.
2. After washing the glass with PBS after fixation, inactivate endogenous peroxidase with methanol containing 0.3% H<sub>2</sub>O<sub>2</sub> at room temperature for 15-30 min. Wash with PBS after inactivation.  
**Note:** 1) The step of peroxidase inactivation is needed only when coloring the section. This process is omitted when performing flow cytometry analysis or only observing with a fluorescence microscope.  
2) When flow cytometric detection is performed subsequently, the above two steps should be done in a conical tube or microtube. The treated cells can be stored for 1-2 months when stored in 70% ethanol at -20°C.
3. Apply 100  $\mu$ l of Permeabilisation Buffer on ice for 2-5 min for well permeation of enzyme reaction mixture. Wash with PBS.
4. Apply 50  $\mu$ l of labeling reaction mixture (consisting of TdT Enzyme 5  $\mu$ l + Labeling Safe Buffer 45  $\mu$ l, prepared and cooled on ice prior to use) on the slide, and incubate in a 37°C humidified chamber for 60-90 min. It is recommended to cover with a coverslip to prevent drying. Terminate the reaction by washing with PBS.  
**Note:** When performing reaction in a tube, mix gently once in every 15 min. to suspend the precipitated cells.

\*For coverslip, Takara Slide Seal for in situ PCR (Cat.#9066, 9067, 9068) is available.

The slides treated as above can be applicable to detection with a fluorescence microscope or flow cytometry. When viewing with a light microscope, follow the procedure described as below.

5. Apply Anti-FITC HRP Conjugate at 37°C for 30 min, and wash with 3-4 times PBS. After coloring with DAB at room temperature for 10-15 min, terminate the reaction by washing with distilled water.  
**Note:** When performing reaction in a tube, sometimes mix gently so that antibody can react uniformly with the cells.
6. Stain the cells with 3% methyl green. Mount and detect with a light microscope.

## [B. Frozen tissue section]

1. Freeze the fresh tissues immediately in an OTC compound. Slice the frozen tissue with a cryostat stick onto a silanized slide. Fix the cells with freshly prepared 4% paraformaldehyde/PBS solution (pH7.4) or acetone at room temperature for 15-30 min. Wash with PBS for 20-30 min.
2. Wash the slides with PBS, and inactivate the endogenous peroxidase using methanol (containing 0.3% H<sub>2</sub>O<sub>2</sub>) at room temperature, 15-30 min.  
**Note:** The step of peroxidase inactivation is needed only when coloring the section. This process is omitted when only observing with a fluorescence microscope.
3. Apply 100  $\mu$  l of Permeabilisation Buffer on ice for 2-5 min. so that enzyme reaction mixture can permeate well.
4. Apply 50  $\mu$  l of labeling reaction mixture (consisting of TdT Enzyme 5  $\mu$  l + Labeling Safe Buffer 45  $\mu$  l, prepared and cooled on ice prior to use) on the slide, and incubate in a 37°C humidified incubator for 60-90 min. It is recommended to cover the slideglass with a coverslip to prevent drying. Terminate the reaction by washing in 3 changes of PBS for 5 min each wash.

The slides treated as above can be applicable to detection with a fluorescence microscope.

5. Apply 70 ml of Anti-FITC HRP Conjugate and incubate at 37°C for 30 min, and wash in 3 change of PBS for 5 min per wash.  
**Note:** 1) Cover the tissue uniformly with the antibody.  
2) It is recommended to cover the slideglass with a coverslip to prevent drying.
6. After coloring with DAB at room temperature for 10-15 min, terminate the reaction by washing with distilled water.
7. Stain with 3% methyl green. Observe using a light microscope after dehydration, penetration and sealing.

\*For coverslip, Takara Slide Seal for in situ PCR (Cat.#9066, #9067, #9068) is available.

## [C. Paraffin embedded tissue section]

1. Deparaffinize the section following the procedure described in "10. Protocol, D. Deparaffinization", page 6. Wash with distilled water. Apply 10~20  $\mu$  g/ml Proteinase K and leave at room temperature for 15 min. Wash with PBS.  
**Note:** when the intensity of staining of apoptosis cells is low, change the reaction condition of Proteinase K (400  $\mu$  g/ml, 5min). When the incubation time is too long, the disruption of the tissue may be caused.
2. Inactivate the endogenous peroxidase by applying 3% H<sub>2</sub>O<sub>2</sub> for 5 min. Wash with PBS.  
**Note:** The step of peroxidase inactivation is needed only when coloring the section. This process is omitted when only observing with a fluorescence microscope.
3. Apply 50  $\mu$  l of labeling reaction mixture (consisting of TdT Enzyme 5  $\mu$  l + Labeling Safe Buffer 45  $\mu$  l, prepared and cooled on ice prior to use) on the slide, and incubate in a 37°C humidified chamber for 60-90 min. It is recommended to cover the slideglass with a coverslip to prevent drying. Terminate the reaction by washing the slides in 3 changes of PBS for 5 min per wash.

The above processes allows detection with a fluorescence microscope.

4. Apply 70  $\mu$  l of Anti-FITC HRP Conjugate and incubate at 37°C for 30 min, and wash in 3 changes of PBS for 5 min. per wash.  
**Note:** 1) Cover the tissue uniformly with the antibody.  
2) It is recommended to cover the slideglass with a coverslip\* to prevent drying.
5. After coloring with DAB at room temperature for 10-15 min, terminate the reaction by washing with distilled water.

6. Stain with 3% methyl green. Observe with a light microscope after dehydration, penetration and sealing.

\* For coverslip, Takara Slide Seal for *in situ* PCR (Cat.#9066, 9067, 9068) is available.

#### [D. Deparaffinization]

1. Apply the followings in order.

Xylene I	for 5 min.
Xylene II	for 5 min.
Xylene III	for 5 min.
100% ethanol	for 5 min.
100% ethanol	for 5 min.
90% ethanol	for 5 min.
80% ethanol	for 5 min.
2. Wash with flowing water for 2 min.
3. Immerse in distilled water.


Please read through this note prior to starting the protocol

#### 11. Note

1. It is recommended to use silanized slideglasses to prevent exfoliation.
2. Humidified chamber should be prewarmed to 37°C .
3. Covering the slides with coverslips during reaction is useful to spread the reaction mixture uniformly (by capillary phenomenon). It also prevents evaporation during incubation.
4. When covering the slides with reaction mixture after washing with PBS, tap off the excess water with filter paper or paper towel.
5. When washing the cells with PBS, pay special attention not to pour PBS directly to the cells because cells exfoliate. When using tissue section as a specimen, wash in 3 changes for 3 min. per wash.

\* For coverslip, Takara Slide Seal for *in situ* PCR (Cat.#9066, 9067, 9068) is available.

## 12. Flow chart of procedure



## 13. Control Slide Experiment

Optical-microscope observation of the stained paraffin-embedded tissue section.



control slide (rat mammary gland)

DAB color development  
Contrast staining 3% methyl green

## 14. Q&A

**Q1:** The intensity of staining of apoptosis cells is low.

**A1:** i) The reaction mixture might not have permeated well in tissues or cells due to steric hindrance. For the enough permeation of the mixture, please adjust the treatment time with Proteinase K or Permeabilization buffer.  
ii) Extend the enzymatic reaction time.  
iii) Extend the antibody reaction time or time of coloring the substrate.

**Q2:** Non-apoptosis cells are stained.

**A2:** Non-specific binding might occur. Repeat the washing steps or add the blocking reagent into the washing buffer, ex. 1%(w/v) BSA, or skimmed milk.

## 15. Related Product

ApopLadder Ex<sup>TM</sup> (Cat.#MK400)

---

**Note:** This product is intended to be used for research purpose only. They are not to be used for drug or diagnostic purposes, nor are they intended for human use. They shall not to be used products as food, cosmetics, or utensils, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please call at +81 77 543 7247 or contact from our website at [www.takara-bio.com](http://www.takara-bio.com).

---