T4 DNA Ligase

Code No. 2011A **Size:** **25,000 U**
Conc.: **350 U/ μ l**

Supplied Reagents:
10X T4 DNA Ligase Buffer **1 ml**

Description:

T4 DNA Ligase catalyzes the formation of phosphodiester bonds between 3'-OH termini and 5'-P termini to join DNA strands. The molecular weight of T4 DNA Ligase is ~62,000 Da. For ligation, the optimal reaction pH is 7.6, and T4 DNA Ligase requires Mg^{2+} and ATP as cofactors. This enzyme can ligate molecules with either cohesive ends or blunt ends.

Storage Buffer:

10 mM	Tris-HCl (pH 7.5)
50 mM	KCl
1 mM	DTT
0.1 mM	EDTA
50%	Glycerol

Storage: -20°C

Source:

Escherichia coli carrying the plasmid that enables highly expression of T4 DNA ligase gene^{2), 3)}

Unit definition:

One unit is the amount of enzyme that ligates more than 90% of 6 μ g of λ DNA-*Hind* III fragments in a 20 μ l mixture in 30 minutes at 16°C. One unit of this enzyme corresponds to 0.008 Weiss units by the ATP-PPi exchange reaction¹⁾ in 66 mM Tris-HCl, pH 7.6 that contains 6.6 mM $MgCl_2$, 10 mM DTT, 66 μ M ATP, and 3.3 μ M [³²P]- $Na_4P_2O_7$.

Reaction mixture for unit definition:

66 mM	Tris-HCl (pH 7.6)
6.6 mM	$MgCl_2$
10 mM	DTT
0.1 mM	ATP
6 μ g/20 μ l	λ DNA- <i>Hind</i> III fragments

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

- 1) Insertion of a DNA fragment into a vector
- 2) Linker or adaptor ligation to a DNA fragment

Note:

The reaction rate of ligation varies depending on the terminal sequence of the DNA molecule. In general, the order is as follows (fastest to slowest):

cohesive-end ligation:

Hind III > *Pst* I > *Eco*RI > *Bam*HI > *Sal*I⁴⁾

blunt-end ligation: *Hae* III > *Alu* I > *Hinc* II > *Sma* I

*Eco*RV > *Sca* I > *Pvu* II > *Nru* I⁴⁾

The ligation rate of *Hind* III end is 10 - 40 times faster than *Sal* I end, and *Hae* III end is 5 - 10 times faster than *Sma* I end.

Composition of Supplied Reagents (store at -20°C):

10X T4 DNA Ligase Buffer	
660 mM	Tris-HCl, pH 7.6
66 mM	$MgCl_2$
100 mM	DTT
1 mM	ATP

References:

- 1) Weiss B, Jacquemin-Sablon A, Live T R, Fareed G C, and Richardson C C. *J Biol Chem.* (1968) **243**: 4543-4555.
- 2) Wilson G G and Murray N E. *J Mol Biol.* (1979) **132**: 471-491.
- 3) Murray N E, Bruce S A, and Murray K. *J Mol Biol.* (1979) **132**: 493-505.
- 4) Hayashi K, Nakazawa M, Ishizaki Y, and Obayashi A. *Nucleic Acids Res.* (1985) **13**: 3261-3271.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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