RNase T1, AOF

Code No. 2152A **Size:** 100 KU
 Conc.: 1,000 U/ μ l

Code No. 2152B **Size:** 500 KU
 Conc.: 1,000 U/ μ l

Description:

RNase T1, AOF (Animal Origin Free) is an endoribonuclease that specifically recognizes and cleaves single-stranded RNA molecules at guanine bases. It hydrolyzes the phosphodiester bond on the 3' side of guanine bases.¹⁾ This enzyme does not require metal ions for its activity and is inhibited by divalent metal ions such as Mg²⁺ and Ca²⁺.²⁾

Product Quality:

1. This product does not contain any human- or animal-derived materials in the final formulation.
2. This product does not contain any β -lactam compounds in the final formulation.

Storage Buffer:

50 mM Tris-HCl, pH 7.5
50% Glycerol

Storage: -20°C

Note: Avoid repeated freeze-thaw cycles.

Source:

Escherichia coli carrying a plasmid containing the gene for Ribonuclease T1 from *Aspergillus oryzae*

Properties:

Molecular mass : approx. 11 kDa
Optimum pH : about pH 7.5

Unit Definition:

One unit is the amount of enzyme that increases the absorbance at 260 nm by 1.0 in 15 minutes at 37°C and pH 7.5, using carrier RNA as the substrate.

Reaction Mixture for Unit Definition:

50 mM Tris-HCl, pH 7.5 (37°C)
2 mM EDTA, pH 7.5
0.2 mg/ml Carrier RNA

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

1. RNA removal from DNA preparations
2. Analysis of Poly(A) length of mRNA

References:

- 1) Takahashi K and Moore S. The Enzymes (ed. Boyer, P.D). *Academic Press, New York*. (1982) Vol. 15B, 3rd ed., 435-468.
- 2) Sato K and Egami F. *J Biochem*. (1957) **44**(11): 753-767.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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- 1) Takahashi K and Moore S. The Enzymes (ed. Boyer, P.D). Academic Press, New York. (1982) Vol. 15B, 3rd ed., 435-468.
- 2) Sato K and Egami F. *J Biochem.* (1957) **44**(11): 753-767.

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