DNA Topoisomerase I

Code No. 2240A **Size:** **100 U**
Conc.: **20 U/ μ l**

Supplied Reagents:

10X DNA Topoisomerase I Buffer **1 ml**
0.1% BSA **1 ml**

Description :

This enzyme has the following four activities³⁾ :

- (1) It relaxes the supercoil of supercoiled circular molecules.
 - (2) It produces knots in single-stranded, circular DNA (knotting), and unwinds knots as well (unknotting).
 - (3) It forms a double-stranded, closed circular DNA from two single-stranded circular DNAs that are complementary to each other.
 - (4) It joins two molecules when a nick exists in one of two double-stranded circular DNAs (catenation), and separates them as well (decatenation).
- The enzyme is obtained from calf thymus, and is active without Mg^{2+} , unlike enzymes from prokaryotes. Prokaryotic DNA topoisomerase I acts only on negatively supercoiled molecules in reaction 1, but this enzyme can relax both positively and negatively supercoiled molecules.

Storage Buffer :

20 mM	Potassium Phosphate (pH 7.2)
50 mM	KCl
0.05 mM	EDTA
1 mM	DTT
0.02 %	BSA
50 %	Glycerol

Storage: -20°C

Source : Calf thymus^{1,2)}

Unit definition :

One unit is the amount of the enzyme that completely relaxes 0.5 μ g of supercoiled pBR322 DNA in 50 μ l of the reaction mixture for 30 minutes at 37 °C.

Reaction mixture for unit definition :

35 mM	Tris-HCl, pH 8.0
72 mM	KCl
5 mM	MgCl ₂
5 mM	DTT
5 mM	spermidine
0.01%	bovine serum albumin
0.5 μ g/50 μ l	supercoiled pBR322 DNA

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Note :

When supercoiled molecules such as pBR322 DNA or ϕ X174 DNA RF I are treated with DNA Topoisomerase and are analyzed with electrophoresis on an agarose gel that contains EtBr, no change in the position of the bands occurs with the enzyme treatment. This is because DNA, which is converted to the relaxed form by DNA Topoisomerase I, resumes its supercoiled configuration when treated with EtBr. Thus, when the activity of DNA Topoisomerase I is to be measured, we recommend that electrophoresis be done with an agarose gel that does not contain EtBr, and then staining with EtBr be done afterwards.

Composition of Supplied Reagents (Store at -20 °C) :

1. 10X DNA Topoisomerase I Buffer
350 mM Tris-HCl (pH 8.0)
720 mM KCl
50 mM MgCl₂
50 mM DTT
50 mM spermidine
2. 0.1% BSA

The supplied buffer has the general composition to be applicable to the reactions to relax the supercoil of supercoiled circular molecules. The supplied buffer gives the same composition as that of used for unit definition. White precipitate may be generated when adding 0.1% BSA directly to the supplied buffer. Prepare the reaction mixture (final conc. 0.01% BSA) by adding the reagents in the order given below. Sterile purified water → Supplied buffer → 0.1% BSA → Substrate DNA

Application example :

pBR322 DNA RFI	0.5 μ g
DNA Topoisomerase I	1 U
10X DNA Topoisomerase I Buffer	2 μ l
0.1% BSA	2 μ l
Sterile purified water	up to 20 μ l

↓
37 °C for 30 min

The supercoiled pBR322 DNA becomes relaxed.

References :

- 1) Ferro A M, Higgins N P, and Olivera B M. *J Biol Chem.* (1983) **258**: 6000-6003.
- 2) Ferro A M and Olivera B M. *J Biol Chem.* (1984) **259**: 547-554.
- 3) Yaginuma K and Koike K. *Protein, Nucleic Acid and Enzyme.* (1983) **28**: 298-305 (in Japanese).

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

DNA Topoisomerase I

Code No. 2240A 容量： 100 U
濃度： 20 U/ μ l

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10 × DNA Topoisomerase I Buffer 1 ml
0.1% BSA 1 ml

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