Alkaline Phosphatase (Calf intestine)

Code No. 2250A **Size:** **1,000 U**
Conc.: **30 U/ μ l**

Supplied Reagents:
10X Alkaline Phosphatase Buffer 1 ml

Description :

This enzyme is a dimer made of subunits each with the molecular weight of 50,000. Its active form is as a metalloenzyme of glycoprotein containing Zn^{2+} . It degrades almost all phosphate monoesters in the same way as the enzyme from *E. coli*, but it does not degrade diphosphate or triphosphate linkages.

Storage buffer :

10 mM	Tris-HCl, pH 8.0
1 mM	MgCl ₂
50 mM	KCl
0.1 mM	ZnCl ₂
50%	Glycerol

Storage: -20°C

Source :

Yeast carrying the plasmid which encodes the gene of calf intestine alkaline phosphatase.

Unit definition :

One unit is the amount of the enzyme that produces 1 μ mol of *p*-nitrophenol per minute at 37°C and pH 9.8, with *p*-nitrophenyl phosphate as the substrate.

Reaction mixture for unit definition :

1 M	diethanolamine, pH 9.8
0.5 mM	MgCl ₂
15 mM	<i>p</i> -nitrophenyl phosphate

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Note :

The optimum pH of this enzyme is near 10 when the concentration of the substrate is high. At low substrate concentrations, the optimum pH is near 8 and the specific activity is less than that at high substrate concentration. This enzyme is quite stable under proper storage conditions. It is irreversibly inactivated (99% or more) after being heated at 65°C for 30 minutes in the presence of a chelating reagent. However, inactivation may not be sufficient, depending on the conditions. For complete inactivation, treatment with phenol is desirable.

Applications :

Dephosphorylation of 5'-terminal phosphate groups from DNA fragments.

Composition of Supplied Reagent :

10X Alkaline Phosphatase Buffer (Store at -20°C)
500 mM Tris-HCl, pH 9.0
10 mM MgCl₂

The above buffer has a general composition applicable to experiments using Alkaline Phosphatase (Calf intestine). It is different from reaction mixture for unit definition.

Application example : Dephosphorylation of DNA

Prepare the reaction mixture by combining the following in a 1.5 ml tube to a total volume of 50 μ l.

DNA fragment in TE Buffer	1 - 20 pmol
Calf Intestine Alkaline Phosphatase (30 U/ μ l)	1 - 2 μ l
10X Alkaline Phosphatase Buffer	5 μ l
Sterile purified water	up to 50 μ l

- Incubate at 37 or 50°C for 30 min
Or incubate at 37°C for 15 min and then incubate at 50°C for 15 min *1
- Phenol/ Chloroform/ Isoamyl alcohol (25 : 24 : 1) extraction (2 times) *2
- Chloroform/ Isoamyl alcohol (24 : 1) extraction
- ← Add 2.5 μ l of 3 M NaCl (final conc. 150 mM)
- Ethanol precipitation (Add 125 μ l (x 2.5 volume) of chilled ethanol, leave at -20°C for 30 - 60 min, and centrifuge.)
- Wash the DNA pellet with 200 μ l of 70% chilled ethanol, and dry it.
- ↓
Dissolve in TE Buffer (< 20 μ l)

- * 1 Single strand DNA and DNAs with 5'-protruding ends can be dephosphorylated at 37°C. But those with blunt or 3'-protruding ends should be incubated at a higher temperature.
- * 2 It is recommended to incubate at 56°C for 30 min in 5 mM EDTA, 0.5% SDS and 50 μ g/ml Proteinase K before phenol extraction to inactivate the enzyme completely. Incubation at 75°C for 10 min prior to phenol extraction is more effective.

References :

- 1) Morton, R K. *Biochem J.* (1957) **65**: 674-682.
- 2) Chappellet-Tordo D, Fosset M, Iwatsubo M, Gache C, and Lazdunski M. *Biochemistry.* (1974) **13**: 1788-1795.
- 3) Besman M, and Coleman J E. *J Biol Chem.* (1985) **260**: 11190-11193.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Alkaline Phosphatase (Calf intestine)

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