# **Recombinant DNase I** (RNase-free)

Code No. 2270A	Size:	1,000 U
	Conc.:	5 U/μl
Supplied Reagent:		
10X DNase I Bu	ıffer	1 ml

## Description:

Recombinant DNase I (RNase-free) is an endodeoxyribonuclease that catalyzes to the same extent the degradation of both single- and double-stranded DNA randomly, and produces 5'-P terminal oligonucleotides. As protease activity and ribonuclease activity have been eliminated, this enzyme is stable around its optimum neutral pH range. So, this enzyme is suitable for RNA preparation at neutral pH.

## Storage Buffer :

Fris-HCl, pH 7.5 (at 4 ℃)
NaCl
CaCl <sub>2</sub>
Glycerol

Storage: -20℃

Recombinant enzyme derived from non-animal host. Source:

## Unit definition :

One unit is the amount of the enzyme that increases the absorbance at 260 nm by 0.001 per minute at 25°C, pH 5.0, with calf thymus DNA as the substrate (Kunitz unit).<sup>1)</sup>

# Reaction mixture for unit definition :

100 mM	Sodium acetate, pH 5.0 (at 25 ℃)
5 mM	MgSO <sub>4</sub>
40 µg/ml	calf thymus DNA

## **Quality Control Data :**

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications:

- 1. For digestion of template DNA without buffer exchange, after *in vitro* transcription with T7 or SP6 RNA Polymerase.<sup>5)</sup>
- 2. For nick translation togather with DNA Polymerase I (Cat.#2130).<sup>2)</sup>
- 3. Making a DNA library for shotgun sequencing in the presence of Mn<sup>2+. 3)</sup>
- 4. For foot-printing analysis of the interaction between DNA and protein.4)
- 5. For digestion of genomic DNA before RT-PCR.

## Note:

As protease activity has been eliminated, this enzyme is stable around its optimum neutral pH range. The enzyme requires a divalent metal ion for its activity. It randomly produces nicks in double-stranded DNA in the presence of Mg<sup>2+</sup>, but in the presence of Mn<sup>2+</sup>, both strands of double stranded DNA are cleaved into fragments. The enzyme loses its activity reversibly with EDTA, and irreversibly by heat treatment at 75°C for 10 min.

#### Composition of Supplied Reagents (Store at -20 °C): 10X DN

DNase I Buf	fer
400 mM	Tris-HCl, pH 7.5
80 mM	MgCl <sub>2</sub>
50 mM	DTT

## Application example :

Digestion of genomic DNA in a sample for RT-PCR

1. Prepare the following reaction mixture

epare are rono migreaction mixturer	
total RNA	20-50 μg
10X DNase I Buffer	5 µ l
Recombinant DNase I (RNase-free)	2 µI(10 U)
Ribonuclease Inhibitor	20 U
DEPC-treated water	up to 50 $\mu$ l
- 1- + - ( +	

2. Incubate for 20 - 30 min at 37°C.

- 3. Perform one of the following procedures to inactivate DNase I.
  - A. Heat treatment
    - (1) Add 2.5 μl of 0.5 M EDTA, and incubate at 80°C for 2 min. (2) Increase reaction volume to 100  $\mu$  l with DEPC treated water.
  - B. Phenol/Chloroform extraction
    - (1) Mix 50  $\mu$  of DEPC treated water and 100  $\mu$  of phenol/ chloroform/isoamyl alcohol (25:24:1) together.
    - (2) Centrifuge at 12,000 rpm for 5 min at room temperature, then transfer the upper layer to a new tube.
    - (3) Add equal amount of chloroform/isoamyl alcohol (24:1) and mix
    - (4) Centrifuge at 12,000 rpm for 5 min at room temperature, then transfer upper layer to new tube.
- 4. Add 10  $\mu$  l of 3 M sodium acetate and 250  $\mu$  l of chilled ethanol, and then mix. Keep it for 20 min at - 80°C.
- 5. Centrifuge at 12,000 rpm for 10 min at 4 °C. Remove the supernatant.
- 6. Wash the precipitate with chilled 70 % ethanol. Centrifuge at 12,000 rpm for 5 min at 4 °C and remove the supernatant.
- 7. Dry the precipitate.
- 8. Dissolve the precipitate in a suitable amount of DEPC-treated water. Confirm the genomic DNA is removed completely by electrophoresis and measure the RNA concentration. When the genomic DNA is not removed completely, increase the amount of enzyme or extend reaction time.

# References:

- 1) Kunitz M. J Gen Physiol. (1950) 33: 349-377.
- 2) Rigby P W J, Dieckmann M, Rhodes C, and Berg P. J Mol Biol. (1977) 113: 237-251.
- 3) Anderson S. Nucleic Acids Res. (1981) 9: 3015-3027. 4) Galas D J and Schmitz A. Nucleic Acids Res. (1978) 5: 3157-3170.
- 5) Green M R, Maniatis T, and Melton D A. Cell. (1983) 32: 681-694.

# Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

# Recombinant DNase I (RNase-free)

Code No. 2270A	容量:	1,000 U
	濃度:	5 U/μl

# 添付試薬:

10 × DNase I Buffer

1 ml

## ● 製品説明

Recombinant DNase I (RNase-free) は、一本鎖および二本鎖の DNA を同 程度にランダムに分解し、5'-P 末端を持つオリゴヌクレオチドを生成さ せる endodeoxyribonuclease である。Mg<sup>2+</sup>存在下では二本鎖 DNA にラ ンダムにニックを入れるが、Mn<sup>2+</sup>存在下では二本鎖の同時切断が起こり、 DNA を断片化する。

本酵素は、実用上問題とならないレベルまで Protease および RNase を除いているため、中性域での安定性が向上しており、RNA の調製の際に安全に用いることができる。

20 mM	Tris-HCl (pH7.5 at 4℃)
50 mM	NaCl
0.1 mM	CaCl <sub>2</sub>
50 %	Glycerol
	50 mM 0.1 mM

●**保存** - 20°C

● 起源 Recombinant enzyme derived from non-animal host.

#### ● 活性の定義

仔牛胸腺 DNA を基質として、25℃、pH5.0 において反応液の 260 nm の 吸光度を 1 分間に 0.001 増加させる酵素活性を 1 U とする (Kunitz unit<sup>1)</sup>)。

## ●活性測定用反応液組成

100 mM	酢酸ナトリウム (pH5.0 at 25℃)
5 mM	MgSO <sub>4</sub>
40 $\mu$ g/ml	基質 DNA

#### ● 品質管理データ

性能試験結果については、各ロットの Certificate of Analysis (CoA) をご 覧ください。CoA はタカラバイオウェブサイトからダウンロードできます。

## ●用途

- 1. T7 あるいは SP6 RNA Polymerase を用いて *in vitro* transcription による RNA 合成を行った後、buffer 置換せずに鋳型 DNA を消化 <sup>5)</sup> する 際に用いられる。
- 2. DNA Polymerase I (製品コード 2130) と共に nick translation<sup>2)</sup> に用 いられる。
- 3. Mn<sup>2+</sup>存在下で shot- gun sequencing のための DNA ライブラリー 作製 <sup>3)</sup> に用いられる。
- 4. DNA- タンパク質の相互作用を調べるフットプリント法<sup>4)</sup>に用いられる。
- 5. RT-PCR 前のゲノム DNA の除去に用いられる。

## ● 使用上の注意

- ・本酵素の至適 pH は中性付近にあり、実用上充分なまで Protease を除 いているために中性付近でも安定である。
- ・本酵素は、活性の発現に二価金属イオンを要求する。
- ・ Mg<sup>2+</sup>存在下では二本鎖 DNA にランダムにニックを入れるが、Mn<sup>2+</sup>存在下では二本鎖の同時切断が起こり、DNA を断片化する。EDTA に より可逆的に失活し、75℃、10 分の熱処理で不可逆的に失活する。

# タカラバイオ株式会社

# ウェブサイト https://www.takara-bio.co.jp

## ● 添付試薬液組成(保存:-20℃)

10 × DNase I Bu	ffer
400 mM	Tris-HCl (pH7.5)
80 mM	MgCl <sub>2</sub>
50 mM	DTT

#### ● 使用例

RT-PCR 前の混入ゲノム DNA 除去のための DNase I 処理方法

1.	以下の反応液を調製する。	
	total RNA	$20 \sim 50 \ \mu g$
	10  imes DNase   Buffer	5 µ l
	Recombinant DNase I (RNase-free)	2 µI (10 Ū)
	RNase Inhibitor	20 U
	DEPC 処理水	50 μl に fill up
h		

- 2. 37℃で 20 ~ 30 分間反応する。
- 以下のいずれかの方法でDNaselを失活させる。
  【熱処理】
  - (1) 2.5µIの0.5 MEDTAを加えて、80℃で2分間インキュベートする。
    (2) DEPC 処理水で 100µIに fill up する。
  - 【フェノール/クロロホルム抽出】
  - (1) 50 µlの DEPC 処理水と 100 µlのフェノール/クロロホルム/ イソアミルアルコール (25:24:1)を加えて混合する。
  - (2) 室温、12,000 rpm で5分間遠心し、上層を新しいチューブに移す。
    (3) 等量のクロロホルム/イソアミルアルコール (24:1) を加えて 混合する。
  - (4) 室温、12,000 rpm で5分間遠心し、上層を新しいチューブに移す。
- 10 µlの3M酢酸ナトリウムと250 µlの冷エタノールを加えて-80℃ で20分間静置する。
- 5. 4℃、12,000 rpm で 10 分間遠心し、上清を捨てる。
- 6. 70%エタノールで沈殿を洗浄し、4 ℃、12,000 rpm で 5 分間遠心し、 上清を捨てる。
- 7. 沈殿を乾燥する。
- 8. 適当量の DEPC 処理水に溶解後、電気泳動によりゲノムDNAが除去できたことを確認し、濃度を測定する。(ゲノム DNA が完全に除去できていない場合は、DNaselの量を増やすか、反応時間を延ばして再度行う。)

## ● 参考文献

- 1) Kunitz M. J Gen Physiol. (1950) 33: 349-377.
- 2) Rigby P W J, Dieckmann M, Rhodes C, and Berg P. *J Mol Biol.* (1977) **113**: 237-251.
- 3) Anderson S. Nucleic Acids Res. (1981) 9: 3015-3027.
- 4) Galas D J and Schmitz A. Nucleic Acids Res. (1978) 5: 3157-3170.
- 5) Green M R, Maniatis T, and Melton D A. Cell. (1983) 32: 681-694.



v202107Da

製品についての技術的なお問い合わせ先 テクニカルサポートライン