Mung Bean Nuclease

Code No. 2420A **Size:** **2,000 U**
Conc.: **40 U/ μ l**

Supplied Reagents:
10X Mung Bean Nuclease Buffer 1 ml

Description:

Mung Bean Nuclease catalyzes the specific degradation of single-stranded regions in single-stranded and double-stranded nucleic acids, and produces mono- and oligonucleotides carrying a 5'-P terminus. An excess of the enzyme is required to degrade double-stranded DNA or RNA and DNA-RNA hybrids, and in this case, AT-rich regions are selectively degraded.²⁾

Storage Buffer:

10 mM Tris-HCl, pH 7.5
0.1 mM zinc acetate
50% glycerol

Storage: -20°C

Source: Mung bean sprouts

Unit definition:

One unit is the amount of the enzyme that converts 1 μ g of heat-denatured calf thymus DNA into acid-soluble form at pH 5.0 in 1 minute at 37°C.

Reaction mixture for unit definition:

30 mM sodium acetate, pH 5.0
100 mM NaCl
1 mM zinc acetate
10% glycerol
0.5 mg/ml substrate DNA

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

1. Mapping of hybridization (S1 mapping). As it has less nibbling activity than S1 Nuclease (Cat. #2410A/B), extra bands are rarely generated.⁵⁾
2. Making the ends of double-stranded DNA blunt-ended (however, the efficiency depends on the base sequence).
3. Cutting out gene coding regions in the presence of formamide.⁶⁾

Note:

This enzyme is stable at pH 7.0. This enzyme loses its activity at pH 5.0 if 0.1 mM Zn²⁺, 1 mM cysteine and 0.005% TritonX-100 is not included in this enzyme. When pH is low, the specificity for single-stranded becomes weak.¹⁾ It becomes completely inactive when heated with EDTA or when treated with 0.01% SDS. Its ability to recognize double-stranded nucleic acids depends on the base sequence. It tends to cleave at ApN and at T(U) pN.³⁾ It completely degrades ApA, but does not degrade G and C. Unlike S1 Nuclease, it does not cleave the strand opposite to that which has been nicked.⁴⁾

Composition of Supplied Reagents:

10X Mung Bean Nuclease Buffer
300 mM sodium acetate, pH 5.0
1,000 mM NaCl
10 mM zinc acetate
50% glycerol

The supplied buffer gives the same composition as that used for the unit definition. The supplied buffer has the general composition to be applicable to the experiments like conversion of protruding termini of double-stranded DNA to blunt ends.

Application Example:

Removal of single-stranded extensions (3' and 5') to leave ligatable blunt ends.

DNA	1 μ g
10X Mung bean Buffer	5 μ l
Mung Bean Nuclease	10 - 30 U
Sterile purified water	up to 50 μ l

↓
— Incubation at 37°C for 10 min
— Phenol/Chloroform extraction
— Ethanol precipitation
— Dissolve in TE Buffer
↓
Use in next reaction

References:

- 1) Kowalski D, Kroeker W D, and Laskowski M, Sr. *Biochemistry*. (1976) **15**: 4457-4467.
- 2) Sheflin L G and Kowalski D. *Nucleic Acids Res.* (1985) **13**: 6137-6154.
- 3) Sung S-C and Laskowski M, Sr. *J Biol Chem.* (1962) **237**: 506-511.
- 4) Kroeker W D and Kowalski D. *Biochemistry*. (1978) **17**: 3236-3243.
- 5) Green M R and Roeder R G. *Cell*. (1980) **22**: 231-242.
- 6) McCutchan T F, Hansen J L, Dame J B, and Mullins J A. *Science*. (1984) **225**: 625-628.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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