# mRNA Cap 2'-O-Methyltransferase

**Code No. 2470A**                      **Size:**                      **2,500 U**  
**Conc.:**                              **50 U/μl**

## Supplied Reagents:

**10X Capping Buffer**                                              **100 μl**  
**S-adenosylmethionine (SAM) (32 mM)**                      **100 μl**

## Description:

mRNA Cap 2'-O-Methyltransferase is an enzyme originating from vaccinia virus and specifically methylates 2'-O position of the first nucleotide of Cap-0 RNA that have a 7-methylguanylate structure in the 5' end of RNA, resulting in the production of Cap-1 RNA. Cap-1 RNA has been known to help to evade innate immune response *in vivo* and promote translation. This product also includes the reaction buffer and S-adenosylmethionine (SAM) as a methyl donor, and allows the efficient production of Cap-1 RNA from Cap-0 RNA. Moreover, the use of Vaccinia Capping Enzyme (Cat. #2460A/B) together with this enzyme allows the production of Cap-1 RNA from uncapped RNA (5'-triphosphate RNA) directly by a single-step reaction.

**Storage:** -20°C

## Source:

*Escherichia coli* carrying a plasmid containing the gene for vaccinia virus mRNA cap 2'-O-methyltransferase

## Properties:

Molecular mass: approx. 40 kDa

## Unit Definition:

One unit is defined as the amount of enzyme that methylates 10 pmol of Cap-0 RNA in 1 hour at 37°C.

## Reaction Mixture for Unit Definition:

1X	Capping Buffer
1 mM	DTT
13.75 μCi/ml	[3H] SAM
2 μg/50 μl	100 nt Cap-0 RNA

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications:

Preparation of Cap-1 mRNA

## Precautions for Use:

1. Don't mix the enzyme vigorously.
2. RNA is less produced or fragmented when reagents, tubes, or pipette tips are contaminated with RNase. Precautions should be taken during experiments to avoid RNase contamination, such as wearing disposable gloves and using tubes and tips exclusively for RNA experiments.

## Protocols:

### 1. Preparation of Cap-1 RNA from Cap-0 RNA

Denatured Cap-0 RNA (10 μg)*1,4	16 μl
10X Capping Buffer	2 μl
SAM (4 mM, dilute 32 mM stock to 4 mM)*2	1 μl
mRNA Cap 2'-O-Methyltransferase (50 U/μl)	1 μl
Total	20 μl*3

Incubate at 37°C for 1 hour. \*4

### 2. Preparation of Cap1-RNA from uncapped RNA (5'-triphosphate RNA)

Denatured RNA transcript (10 μg)*1,4	14 μl
10X Capping Buffer	2 μl
GTP (10 mM)	1 μl
SAM (4 mM, dilute 32 mM stock to 4 mM)*2	1 μl
Vaccinia Capping Enzyme (10 U/μl)	1 μl
mRNA Cap 2'-O-Methyltransferase (50 U/μl)	1 μl
Total	20 μl*3

Incubate at 37°C for 1 hour. \*4

\*1 In order to remove secondary structure on the 5' end of the transcript, heat-denature the RNA prior to the reaction.

(1) Take 10 μg of the RNA and adjust the volume to 16 μl (Cap-1 RNA prep from Cap-0 RNA) or 14 μl (Cap-1 RNA prep from uncapped RNA) by RNase-free water.

(2) Heat-denature at 65°C for 5 - 10 min, immediately followed by an incubation on ice for 5 min.

\*2 SAM is unstable. Dilute necessary amount of 32 mM stock solution by RNase-free water prior to the reaction, and keep it on ice until use.

\*3 It can be scaled up depending on the amount required.

\*4 If the RNA length is less than 200 nt, extend the incubation time to enhance the methylation efficiency (and capping efficiency) (1 hr → 2 hrs).

Note: RNase inhibitor can be used to enhance the stability of RNA during the reaction. Add at a final concentration of 1 U/μl if necessary.

## References:

- 1) Barbosa E and Moss B. *J Biol Chem.* (1978) **253**: 7698-7702.
- 2) Kuge H, Brownlee G G, Gershon P D, and Richter J D. *Nucleic Acids Res.* (1998) **26**: 3208-3214.
- 3) Hyde J L and Diamond M S. *Virology.* (2015) **479-480**: 66-74.

## Related products:

Takara IVTpro™ mRNA Synthesis System (Cat. #6141)

Cloning Kit for mRNA Template (Cat. #6143)

Takara IVTpro™ T7 mRNA Synthesis Kit (Cat. #6144)

T7 RNA Polymerase ver.2.0 (Cat. #2541A)

Vaccinia Capping Enzyme (Cat. #2460A/B)

Recombinant RNase Inhibitor (Cat. #2313A/B)

RNase-free Water (Cat. #9012)

NucleoSpin RNA Clean-up (Cat. #740948.10/50/250)\*

\* Not available in all geographic locations. Check for availability in your area.

IVTpro is a trademark of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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Denatured Cap-0 RNA (10 μg)*1,4	16 μl
10×Capping Buffer	2 μl
SAM (4 mM, dilute 32 mM stock to 4 mM)*2	1 μl
mRNA Cap 2'-O-Methyltransferase (50 U/μl)	1 μl
<hr/>	
Total	20 μl*3

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### 2. Uncapped RNA (5'-triphosphate RNA) から Cap-1 RNA の調製

Denatured RNA transcript (10 μg)*1,4	14 μl
10×Capping Buffer	2 μl
GTP (10 mM)	1 μl
SAM (4 mM, dilute 32 mM stock to 4 mM)*2	1 μl
Vaccinia Capping Enzyme (10 U/μl)	1 μl
mRNA Cap 2'-O-Methyltransferase (50 U/μl)	1 μl
<hr/>	
Total	20 μl*3

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