Wide-Range DNA Ladder (50-10,000 bp) Ver.2

Code No. 3418A **Size:** 500 μ l
(for 100 lanes)
Conc.: 680 ng/5 μ l

Description:

This product is a DNA molecular size marker consisting of 16 fragments ranging from 50 bp to 10,000 bp. To facilitate easy identification of fragment sizes after gel electrophoresis, the 500, 1,000, 2,000, and 4,000 bp fragments are designed to appear with higher intensity when stained with EtBr or SYBR™ Green I Nucleic Acid Gel Stain (Cat. #5760A/5761A)*. The product is premixed with loading dyes and glycerol and is ready for direct loading onto an agarose gel. The recommended loading volume is 5 μ l per lane. Each 5 μ l aliquot contains approximately the following amounts of DNA.

| Fragment size (bp) | DNA amount (ng/5 μ l) |
|--------------------|---------------------------|
| 10,000 | 30 |
| 8,000 | 30 |
| 6,000 | 30 |
| 4,000 | 80 |
| 3,000 | 30 |
| 2,000 | 80 |
| 1,600 | 30 |
| 1,400 | 30 |
| 1,000 | 80 |
| 750 | 30 |
| 500 | 80 |
| 400 | 30 |
| 300 | 30 |
| 200 | 30 |
| 100 | 30 |
| 50 | 30 |

Form:

| | |
|---------|------------------|
| 8.33 mM | Tris-HCl, pH 8.0 |
| 5.83 mM | EDTA |
| 6% | Glycerol |
| 0.0083% | Bromophenol Blue |
| 0.0058% | Xylene Cyanol |

Storage:

-20°C

Note: After thawing, store at 4°C and use as soon as you can to avoid repeated freeze-thaw cycles and minimize contamination risk.

Applications:

Used as a DNA molecular size marker in agarose gel electrophoresis.

Application Example:

1. Load 5 μ l of this product onto an agarose gel and run electrophoresis, for example using 1% Agarose L03 "TAKARA" (Cat. #5003/5003B)*.
2. Stain the gel with EtBr or SYBR Green I Nucleic Acid Gel Stain (Cat. #5760A/5761A)*.

* Not available in all geographic locations. Check for availability in your area.

SYBR is a trademark of Life Technologies Corporation.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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