# pAUR123 DNA

**Code No. 3602**

**Size: 20  $\mu$ g**

**Conc.: 1  $\mu$ g/ $\mu$ l**

\* 2 years from date of receipt under proper storage conditions.


**Description:**

pAUR vectors have a mutant *AUR1-C* gene derived from *S. cerevisiae* as a selective marker that confers Aureobasidin A (AbA)-resistance on cells for yeast transformation system. pAUR123 is a protein expression vector and *E. coli*-yeast shuttle vector derived from pAUR112. It contains a constitutive promoter of ADH1 gene<sup>3</sup> (Alcohol dehydrogenase 1 gene) for protein expression in yeast cells.  
(Note: This vector is not suitable for high expression of protein.)

**Form:** 10 mM Tris-HCl, pH 8.0  
1 mM EDTA

**Storage:** -20°C

**Vector map of pAUR123 DNA:**

- AUR1-C*: AbA-resistant gene in *S. cerevisiae*
- ARS: Replication origin in *S. cerevisiae*
- CEN: Centromere in *S. cerevisiae*
- Amp<sup>r</sup>: Selective marker in *E. coli*
- Ori: Replication origin in *E. coli*
- PADH1: Promoter of ADH1 gene
- TADH1: Terminator of ADH1 gene

**GenBank:** Accession No. AB012284

**Preparation:** cccDNA was purified by CsCl-EtBr ultracentrifugation

**Base pairs:** 6,982 bp

**Quality Control Data:**

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.


**Usage:**

This is a vector for transformation of yeast *S. cerevisiae* using antibiotic Aureobasidin A. (shuttle vector for protein expression)

**References:**

- 1) Takesako K, Kuroda H, Inoue T, Haruna F, Yoshikawa Y, Kato I, Uchida K, Hiratani T, and Yamaguchi H. *J Antibiot.* (1993) **46**: 1414-1420.
- 2) Hashida-Okado T, Ogawa A, Endo M, Yasumoto R, Takesako K, and Kato I. *Mol Gen Genet.* (1996) **251**: 236-244.
- 3) Mumberg D, Muller R, and Funk M. *Gene.* (1995) **156**: 119-122.

**Cloning site of pAUR123:** bold: stop codon

**Note**

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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