# pAUR316 DNA

Code No. 3605

Size: 20 µg

Conc.: 0.5 µg/µl

\* 2 years from date of receipt under proper storage conditions.

## Description:

pAUR316 is an *E. coli* - *Aspergillus* shuttle vector which can replicate autonomously in *A. nidulans*. This vector has Aureobasidin A (AbA) resistant gene, *aurA*<sup>r</sup>, derived from *A. nidulans* as a selective marker, AMA1 region of *A. nidulans* as a replication origin in *Aspergillus* cells, and ampicillin resistant gene, *Amp*<sup>r</sup>, as a selective marker in *E. coli*. The restriction endonuclease sites of *Bam*HI, *Aor*51HI, *Afl*II, *Bgl*II, and *Sse*8387I can be used as a cloning site.


**Form:** 10 mM Tris-HCl, pH 8.0  
1 mM EDTA

**Storage:** -20°C

**Preparation:** cccDNA was purified by CsCl-EtBr ultracentrifugation.

**Base pairs:** 11,663 bp

## Vector map of pAUR316:

*aurA*<sup>r</sup>: AbA - resistant gene in *A. nidulans*  
AMA1: Replication oregon in *A. nidulans*  
*Amp*<sup>r</sup>: Selective marker in *E. coli*  
Ori: Replication origin in *E. coli*

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.


## Application:

This is a vector for transformation using antibiotic AbA of *Aspergillus* species which has the AbA-sensitivity, such as *A. nidulans*. It is maintained as a plasmid in *Aspergillus* cells.

## References:

- 1) Takesako K, Kuroda H, Inoue T, Haruna F, Yoshikawa Y, Kato I, Uchida K, Hiratani T, and Yamaguchi H. *J Antibiot.* (1993) **46**: 1414-1420.
- 2) Kuroda M, Hashida-Okado T, Yasumoto R, Gomi K, Kato I, and Takesako K. *Mol Gen Genet.* (1999) **261**: 290-296.

## Cloning site of pAUR316:

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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