Pfu Pyroglutamate Aminopeptidase

Code No. 7334

Size: 10 mU (at 37°C)

Supplied Reagent:
Pfu-PGAP 5X Buffer

1 ml

Description :

Pfu Pyroglutamate Aminopeptidase liberates the N-terminal pyroglutamic acid from proteins and peptides.

Source :

Escherichia coli carrying plasmids encoding the *Pyrococcus furiosus* pyroglutamate aminopeptidase gene

Form :

Lyophilized
Reconstitute with 50 μ l of diluted supplied buffer (1X).

Storage :

-20°C
Reconstituted solution should be stored at 4°C.
For long-term storage, it should be stored in aliquots at -20°C.
Avoid repeated freeze-thaw cycles.

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Definition of Activity :

One unit of enzyme activity corresponds to the amount required to hydrolyze 1 μ mol of Pyroglutamate-*p*-nitroanilide at pH 7.0, in one minute, at 37°C.

Supplied Buffer (5X) :

Pfu-PGAP 5X Buffer ;

250 mM Sodium phosphate (pH 7.0)
50 mM DTT
5 mM EDTA

Properties :

Molecular weight : 24,072 (amino acid composition)
Approx. 28,000 (SDS-PAGE)
Optimum temperature : 95 - 100°C
Thermostability : Approximately 90% activity is confirmed when reacted at pH 7.0 for 150 min at 75°C.
Optimum pH : pH 6.0 - 9.0
Stable pH range : pH 5.0 - 9.0
 \geq 80% activity is confirmed in this range when reacted for 30 min at 75°C.
Tolerance to denaturants :
Stable against \leq 1 M Urea
 \leq 1 M guanidine-HCl
 \leq 0.01% SDS
Stabilizers : DTT, EDTA
Inhibitors : PCMB, Hg²⁺

References :

- 1) Hamazume Y, Mega T, and Ikenaka T. *J Biochem.* (1987) **101**: 217-223.
- 2) Shimada Y, Sugihara A, Tominaga Y, Izumi T, and Tsunasawa S. *J Biochem.* (1989) **106**: 383.
- 3) Tsunasawa S, Nakura S, Tanigawa T, and Kato I. *J Biochem.* (1998) **124**: 778-783.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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