

For Research Use

TAKARA

***E. coli* JM109 Electro-Cells**

Product Manual

Table of Contents

I.	Description.....	3
II.	Components	3
III.	Storage	3
IV.	Protocol.....	4
V.	Quality.....	6
VI.	Genotype.....	6
VII.	Cell Density	6
VIII.	References.....	6
IX.	Related Product	6

I. Description

Electro-Cells are specially prepared by Takara Bio to be best appropriate for electroporation method. Electroporation method is used to transfer DNA into a cell by breaking cytoplasmic membrane by high voltage pulse. As this electro-cells offers high transformation efficiency and good reproducibility, it is especially useful in transferring small amount of sample into a *E. coli* in less time.

As *E. coli* JM109 Electro-Cells contains F' plasmid, it can be used as a host of M13 phage vector as well as for preparation of DNA library or subcloning. When transformation of pUC vectors or transduction of M13 phage vector DNAs, recombinants can be selected easily by X-Gal and IPTG utilizing the α -complimentarity to β -galactosidase of the electro-cells.

X-Gal : 5-Bromo-4-Chloro-3-Indolyl- β -D-Galactoside

IPTG : Isopropyl- β -D-thiogalactopyranoside

II. Components

<i>E. coli</i> JM109 Electro-Cells	50 μ l x 10
pUC19 plasmid (10 pg/ μ l)	10 μ l
SOC Medium*	1 ml x 10

* SOC Medium:	2%	Tryptone
	0.5%	Yeast extract
	10 mM	NaCl
	2.5 mM	KCl
	10 mM	MgSO ₄
	10 mM	MgCl ₂
	20 mM	Glucose

III. Storage

-80°C

Note : Store at -80°C or lower. If the storage temperature is not maintained consistently, the transformation efficiency will be reduced. You may determine the transformation efficiency of stored cells by using the included pUC19 control.

Do not store in liquid nitrogen.

IV. Protocol

A. Transformation of a plasmid vector

- (1) Thaw *E. coli* JM109 Electro-Cells (50 μ l) in an ice bath just before use.
- (2) Add 1-2 μ l of DNA solution*¹ into the thawed cell suspension.
- (3) Transfer the mixture of cells and DNA to a cold 0.1 cm electroporation cuvette.
- (4) After applying pulse*², immediately add 1 ml of SOC Medium (precooled in an ice bath).
- (5) Incubate by shaking (160 - 225 rpm) for 1 hour at 37°C.
- (6) Plate on selective media. Less than 100 μ l is applied to a ϕ 9 cm plate.
- (7) Incubate overnight at 37°C.
 - * 1 When sample DNA solution contains salt, dilute with TE buffer or sterile purified water. Or desalting by ethanol precipitation is recommended (< 10 ng is recommended).
 - * 2 Takara Bio uses BIO-RAD MicroPulser, and the electrical condition is 1.5 kV. In the case of BIO-RAD Gene Pulser, standard electrical conditions are 200 Ω , 25 μ F and 1.5 kV.

B. Transduction of a M13 phage vector

- (1) Follow the step (1)-(4) mentioned in IV-A.
- (2) Add 200 μ l of the host (*E. coli* JM109, A_{600} =0.8 - 1.0) into 3 ml of YT soft agar (preincubated at 46 - 48°C).
- (3) Add a proper amount of the solution prepared at (1) into the agar, mix, and immediately spread it onto a L-plate.
- (4) Incubate at room temperature for 10 - 15 min and then at 37°C over night.

[Please read before proceeding]

1. Place a vial of electro-cells in a dry ice / EtOH bath immediately upon removal from -80°C freezer. Keep cells in bath until you are ready to proceed.
2. When using 50 μ l of electro-cells, apply high-purified sample DNA in less than 10 ng. If not, transformation efficiency might decrease.
3. If using large size of DNA (> 7 kb) transformation efficiency might decrease.
4. When changing an experiment scale, optimum condition should be considered.
5. L-broth or ϕ b-broth can be used instead of SOC Medium. In this case, lower efficiency might be obtained.

<u>L-broth</u> :	<u>Ingredient</u>	<u>per liter water</u>
	Tryptone	10 g
	Yeast extract	5 g
	NaCl	5 g

Adjust to around pH 7.5 with 1 N NaOH and autoclave.

<u>ϕ b-broth</u> :	<u>Ingredient</u>	<u>per liter water</u>
	Tryptone	20 g
	Yeast extract	5 g
	MgSO ₄ ·7H ₂ O	5 g

Adjust to around pH 7.5 with 1 N KOH and autoclave.

6. When diluting, use SOC Medium which has been added in the step (4) of A.

7. <u>YT soft agar</u> :	<u>Ingredient</u>	<u>per 100 ml water</u>
	Tryptone	0.8 g
	Yeast extract	0.5 g
	NaCl	0.5 g

Adjust to around pH 7.6 with 1 N NaOH, add agar to the concentration of 0.6%, and autoclave.

8. Host can be prepared by culturing electro-cells.
When DNA is inserted in a cloning site of M13 phage vector, recombinants can be easily selected by adding X-Gal and IPTG into YT-soft agar because un-recombinants appears in blue plaque.

9. <u>YT-plate</u> :	<u>Ingredient</u>	<u>per liter water</u>
	Tryptone	8 g
	Yeast extract	5 g
	NaCl	5 g

Adjust to around pH 7.5 with 1 N NaOH, add agar to the concentration of 1.5%, and autoclave.

10. When adding X-Gal or IPTG, follow the procedures described as below:
 - Add 100 mM IPTG at 100 - 300 μ l/100 ml agar medium and 25 μ l/3 ml soft agar.
 - Add 20 mg /ml X-Gal (dissolved in dimethylformamide) at a ratio of 200 - 300 μ l/100 ml agar medium, and 50 μ l/3 ml soft agar.
11. It is not recommended to freeze and store the thawed electro-cells. However, if necessary, freeze in a dry ice/EtOH bath and return to -80°C. The transformation efficiency can be lowered by more than one magnitude.

V. Quality

10 pg of pUC19 was transformed and selected by Amp⁺ selective media plating.
Transformation efficiency > 1 x 10⁹ cfu/μg pUC19

VI. Genotype

E. coli JM109 : *recA1, endA1, gyrA96, thi-1, hsdR17*(r_k-m_k⁺), e14⁻(*mcrA*⁻)
supE44, relA1, Δ (lac-proAB)/F⁺[*traD36, proAB*⁺, *lacI*^q, *lacZ* Δ M15]

VII. Cell density

1 x 10¹⁰ bacteria/ml

VIII. References

- 1) Dower W J, Miller J F, and Ragsdale C W. *Nucl Acids Res.* (1988) **16**: 6127.
- 2) Bottger E C. *Biotechniques.* (1988) **6**: 878.

IX. Related product

E. coli JM109 Competent Cells (Cat. #9052)

NOTE: This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.
