For Research Use

TakaRa

E. coli HB101 Competent Cell

Product Manual

v201807Da

Cat. #9051 v201807Da

Table of Contents

Components	3
, ,	
	Components Storage Description Protocols Transformation Efficiency Genotype Cell Density Reference



Cat. #9051 v201807Da

I. Components

<i>E. coli</i> HB101 Competent Cells pBR322 plasmid (0.1 ng/ µ l) SOC Medium *		100 μl x 10 10 μl 1 ml x 10
* SOC Medium :	2 % 0.5 % 10 mM 2.5 mM 10 mM 10 mM 20 mM	Tryptone Yeast extract NaCl KCl MgSO4 MgCl ₂ Glucose

II. Storage

-80°C

Note: If it is not stored at -80°C, transformation efficiency may decrease. In this case, we recommend confirming the efficiency with the supplied pBR322 plasmid prior to using it in an application. Do not store in liquid nitrogen.

III. Description

Competent cells are prepared by Hanahan's method modified by Takara Bio and have a transformation efficiency of >1 x 10⁸ cfu/ μ g when 100 μ l of the cells are transformed by 1 ng of pBR322 plasmid.

E. coli HB101 Competent Cells can be used for preparation of a DNA library or subcloning of a recombinant plasmid.

IV. Protocols (Transformation of a Plasmid Vector)

- 1) Thaw *E. coli* HB101 Competent Cells in an ice bath just before use.
- 2) After thawing, mix the competent cells gently, and transfer 100 μ l of the cells into a 14-ml, round-bottom tube (e.g., Falcon tube, etc.).
- 3) Add < 10 ul of plasmid DNA (\leq 10 ng is recommended).
- 4) Keep the tube in the ice bath for 30 min.
- 5) Incubate the tube for 45 sec at 42°C.
- 6) Return the tube to the ice bath for 1 2 min.
- 7) Add SOC Medium (pre-incubated at 37° C) up to a final volume of 1 ml.
- 8) Incubate by shaking (160 225 rpm) for 1 hour at 37°C.
- 9) Plate the cell suspension on selective media.*
- 10) Incubate overnight at 37°C.
 - * Plate no more than 100 μ l for a 9-cm plate.



Takara

[Precautions for use]

- 1) Transfer a vial of competent cells immediately to a dry ice / EtOH bath from a -80°C freezer. Keep the cells in the bath until you are ready to proceed.
- 2) For transformation, you may use 1.5-ml microcentrifuge tubes instead of 14-ml, round-bottom tubes (BD Code: 352059 or 352057, etc.), but the transformation efficiency may be reduced.
- 3) When using 100 μ l of competent cells, use less than 10 ng of highly purified plasmid DNA. Otherwise, transformation efficiency might decrease.
- 4) When changing an experiment's scale, optimum conditions should be considered.
- 5) L-broth or φ b-broth can be used instead of SOC Medium. However, transformation efficiency might be slightly decreased.

• <u>L-broth :</u>	Ingredient	Per liter water		
	Tryptone	10 g		
	Yeast extract	5 g		
	NaCl	5 g		
Adjust to \sim pH 7.5 with 1N NaOH and autoclave.				
$\cdot \varphi$ b-broth:	Ingredient	Per liter water		
	Tryptone	20 g		
	Yeast extract	5 g		
	MgSO4·7H ₂ O	5 g		

Adjust to ~ pH 7.5 with 1N KOH and autoclave.

- 6) When diluting the transformation mixture, use the medium added in the Step 7 of Section IV.
- 7) Once the competent cells have been thawed, refreezing for storage is not recommended. If this is unavoidable, freeze the cells quickly on dry ice/ethanol and store them promptly at -80°C. However, the transformation efficiency will be lowered by at least one order of magnitude.

V. Transformation Efficiency

1 ng of pBR322 plasmid was transformed and the transformants were selected on an Amp+ selective medium plate. Transformation efficiency \geq 1 x 10⁸ cfu / µg pBR322 plasmid

VI. Genotype

E. coli HB101: F⁻, *hsd* S20(rB-, mB-), *recA*13, *ara*-14, *proA*2, *lacY*1, *galK*2, *rpsL*20 (str), *xyl*-5, *mtl*-1, *supE*44, *leuB*6, *thi*-1.

VII. Cell Density

1 - 2 x 10⁹ bacteria/ml

VIII. Reference

Hanahan D. J Mol Biol. (1983)166: 557.

NOTE : This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.



Cat. #9051

v201807Da