

Cat. # 9055

For Research Use

TAKARA

***E. coli* MV1184 Competent Cells**

Product Manual

v202010Da

Table of Contents

| | | |
|-------|-----------------------|---|
| I. | Description | 3 |
| II. | Components..... | 3 |
| III. | Storage..... | 3 |
| IV. | Protocol | 4 |
| V. | Quality..... | 6 |
| VI. | Genotype..... | 6 |
| VII. | Cell Density | 6 |
| VIII. | Reference | 6 |
| IX. | Related Products..... | 6 |

I. Description

E. coli MV1184 Competent Cells are prepared by Hanahan's method modified by Takara Bio and have a transformation efficiency of $>1 \times 10^7$ cfu/ μ g when 100 μ l of the cells are transformed with 1 ng of pUC119. As *E. coli* MV1184 is an amber suppressor(-) strain in which only non-amber DNA vectors can propagate, this product can be used for selecting amber-mutated DNA. As containing F' plasmid, it can be also used as a host of M13 phage vectors or phagemid vectors in preparation of ssDNA. When transformation of pUC vectors or transduction of M13 phage vector DNAs, recombinants can be selected easily by adding X-Gal and IPTG to a media utilizing α -complementarity to β -galactosidase of the Competent Cells.

X-Gal : 5-Bromo-4-Chloro-3-Indolyl- β -D-Galactoside

IPTG : Isopropyl- β -D-thiogalactopyranoside

II. Components

| | |
|---------------------------------------|------------------|
| <i>E. coli</i> MV1184 Competent Cells | 100 μ l x 10 |
| pUC119 plasmid (0.1 ng/ μ l) | 10 μ l |
| SOC Medium* | 1 ml x 10 |

| | |
|--------------|-------------------------|
| * SOC Medium | 2% Tryptone |
| | 0.5% Yeast extract |
| | 10 mM NaCl |
| | 2.5 mM KCl |
| | 10 mM MgSO ₄ |
| | 10 mM MgCl ₂ |
| | 20 mM Glucose |

III. Storage

-80°C

Note: If it is not stored at -80°C, transformation efficiency may decrease. In this case, it is recommended to confirm the efficiency by using supplied pUC119 prior to use an application.

Do not store in liquid nitrogen.

IV. Protocol

A. Transformation with a plasmid vector

- (1) Thaw *E. coli* MV1184 Competent Cells on ice bath just before use.
- (2) Gently mix cells and transfer 100 μ l of Competent Cells into a 14 ml round-bottom tube.

Important : Do not use a vortex to mix cells.

- (3) Add DNA sample (\leq 10 ng is recommended).
- (4) Keep in the ice bath for 30 min.
- (5) Incubate cells for 45 sec at 42°C.
- (6) Return to the ice bath for 1 - 2 min.
- (7) Add SOC Medium (pre-incubated at 37°C) up to a final volume of 1 ml.
- (8) Incubate by shaking (160 - 225 rpm) for 1 hour at 37°C.
- (9) Plate on selective media*.
* 100 μ l or less is recommended for plating on ϕ 9 cm dish.
- (10) Incubate overnight at 37°C.

B. Transduction with a M13 phage vector

- (1) Follow step (1) - (8) of IV-A.
- (2) Add 200 μ l of the host cells (*E. coli* MV1184, A_{600} =0.8 - 1.0) into 3 ml of YT soft agar (pre-incubated at 46 - 48°C).
- (3) Add a proper amount of the solution prepared at step (1) into the agar, mix, and immediately spread it onto a YT-plate.
- (4) Incubate at room temperature for 10 - 15 min and then, at 37°C overnight.

[Please read before proceeding]

1. Place a vial of the Competent Cells in a dry ice / EtOH bath immediately upon removal from -80°C freezer. Keep cells in bath until you are ready to proceed.
2. You may use 1.5 ml microcentrifuge tubes instead of 14 ml round-bottom tubes (CORNING Code: 352059 or 352057, etc.) for transformation, but it may reduce efficiency.
3. When using 100 μ l of Competent Cells, apply 10 ng or less of highly purified plasmid DNA. If not, transformation efficiency might decrease.
4. When changing an experiment scale, optimum condition should be considered.
5. L-broth or ϕ b-broth can be used instead of SOC Medium. In this case, lower efficiency might be obtained.

| <u>L-broth:</u> | <u>Ingredient</u> | <u>per 1 L water</u> |
|-----------------|-------------------|----------------------|
| | Tryptone | 10 g |
| | Yeast extract | 5 g |
| | NaCl | 5 g |

Adjust to around pH 7.5 with 1 N NaOH and autoclave.

| <u>ϕ b-broth:</u> | <u>Ingredient</u> | <u>per 1 L water</u> |
|-----------------------------------|---------------------------------------|----------------------|
| | Tryptone | 20 g |
| | Yeast extract | 5 g |
| | MgSO ₄ · 7H ₂ O | 5 g |

Adjust to around pH 7.5 with 1 N KOH and autoclave.

6. Medium contents

| <u>YT soft agar:</u> | <u>Ingredient</u> | <u>per 100 ml water</u> |
|----------------------|-------------------|-------------------------|
| | Tryptone | 0.8 g |
| | Yeast extract | 0.5 g |
| | NaCl | 0.5 g |

Adjust to around pH 7.6 with 1 N NaOH, add agar to the concentration of 0.6%, and autoclave.

| <u>YT-plate:</u> | <u>Ingredient</u> | <u>per liter</u> |
|------------------|-------------------|------------------|
| | Tryptone | 8 g |
| | Yeast extract | 5 g |
| | NaCl | 5 g |

Adjust to pH 7.6 with 1 N NaOH, add agar to the concentration of 1.5%, and autoclave.

7. *E. coli* MV1184 cells can be prepared by culturing Competent Cells.
8. When adding X-Gal or IPTG, follow the procedures described as below:
 - Add 100 mM IPTG to be 100 μ l/100 ml agar medium and 25 μ l/3 ml soft agar.
 - Add 20 mg X-Gal/ml dimethylhormeamide to be 200 μ l/100 ml medium and 50 μ l/3 ml soft agar.
9. It is not recommended to refreeze and store the thawed Competent Cells. However, if necessary, freeze in a dry ice/EtOH bath and return to -80°C. The transformation efficiency can be lowered by more than one magnitude.

V. Quality

1. Transformation efficiency
When 100 μ l of *E. coli* MV1184 Competent Cells were transformed with 1 ng of pUC119 according to IV-A and the transformants were selected on a LB-agar plate containing Ampicillin, the transformation efficiency was $> 1.0 \times 10^7$ cfu/ μ g pUC119.
2. Stability of F' plasmid
The white colonies were less than 1% of the total colonies when *E. coli* MV1184 Competent Cells were transformed with pUC119 and plated on a LB-agar plate containing Ampicillin (100 μ g/ml), IPTG (0.3 mM), and X-Gal (60 μ g/ml).

VI. Genotype

E. coli MV1184 : *ara*, Δ (*lac-proAB*), *rpsL*, *thi* (Φ 80 *lacZ* Δ M15), Δ (*srl-recA*) 306 : *Tn10* (*tet^r*)/F' [*traD36*, *proAB⁺*, *lac^R*, *lacZ* Δ M15]

VII. Cell density

1 - 2 x 10⁹ bacteria/ml

VIII. Reference

Hanahan D. *J Mol Biol.* (1983) **166**: 557.

IX. Related products

pUC118 DNA (Cat. #3318)
pUC119 DNA (Cat. #3319)
pUC118 DNA/BAP (Cat. #3320 - 3324)
pTV118N DNA (Cat. #3328)
X-Gal (5-Bromo-4-Chloro-3-Indolyl- β -D-Galactoside) (Cat. #9031)
IPTG (Isopropyl- β -D-thiogalactopyranoside) (Cat. #9030)

NOTE: This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.
