# Recombinant RNase Inhibitor GMP grade

Code No. T2431S001	Size: 40 KU Conc.: 40 U/ $\mu$ l
Code No. T2431S010	Size: 400 KU Conc.: 40 U/ $\mu$ l
Code No. T2431S100	Size: 4,000 KU Conc.: 40 U/ $\mu$ l

## Description:

Recombinant RNase Inhibitor GMP grade is a variant of the ribonuclease inhibitor derived from porcine liver. While this product maintains characteristics similar to RNase inhibitors from porcine liver and human placenta,<sup>1,2</sup> the stability of this protein has been increased by introducing a mutation in one of the cysteine residues, which is responsible for the denaturation due to its easy oxidation.<sup>3,4</sup> This product can be directly added in various reaction mixtures such as *in vitro* transcription and capping reaction, where the integrity of RNA is indispensable.

## Quality Statement:

1. This product is manufactured in accordance with relevant GMP guidelines. For more information, please refer to our website.
2. No animal- or human-derived components are used in the manufacture of this product.
3. This product uses primary materials that do not contain  $\beta$ -lactam compounds.

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Storage Buffer:

20 mM	HEPES-KOH, pH 7.5
50 mM	KCl
5 mM	DTT
50%	Glycerol

**Storage:** -20  $\pm$  5°C

## Source:

*Escherichia coli* carrying a plasmid containing the gene for a variant of ribonuclease inhibitor from porcine liver

## Properties:

Molecular mass : approx. 52 kDa  
Optimum pH : Maximum at pH 7 - 8, though the inhibition is observed in a broad pH range

## Unit Definition:

One unit is defined as the amount of protein that inhibits 50% of RNase A activity, when fluorescent RNA probes are used as the substrate for 5 ng of RNase A.

## Applications:

1. *In vitro* transcription/translation (1 U/ $\mu$ l reaction)<sup>6,7)</sup>
2. Capping reaction (1 U/ $\mu$ l reaction)

\* The parentheses show the examples of RNase inhibitor concentration in each reaction mixture.

## Precautions for Use:

1. Don't mix the protein vigorously.
2. Include DTT in the reaction to maintain the activity, as far as its concentration doesn't interfere with the reaction (e.g., 1 mM).

## References:

- 1) Burton L E and Fucci N P. *Int J Pept Protein Res.* (1982) **19**: 372-379.
- 2) Blackburn P, Wilson G, and Moore S. *J Biol Chem.* (1977) **252**: 5904-5910.
- 3) Kim B-M, Schultz L W, and Raines R T. *Protein Sci.* (1999) **8**: 430-434.
- 4) Dickson K A, Haigis M C, and Raines R T. *Prog Nucleic Acid Res Mol Biol.* (2005) **80**: 349-374.
- 5) Blackburn P. *J Biol Chem.* (1979) **254**: 12484-12487.
- 6) Eichler D C, Tatar T F, and Lasater L S. *Biochem Biophys Res Commun.* (1981) **101**: 396-403.
- 7) Scheele G and Blackburn P. *Proc Natl Acad Sci USA.* (1979) **76**: 4898-4902.

## Related Products:

T7 RNA Polymerase GMP grade (Cat. #T2411S001/S010/S100)  
PrimeCap™ T7 RNA Polymerase (low dsRNA) GMP grade  
(Cat. #T2422S001/S010/S100)  
10X IVT Reaction Buffer GMP grade (Cat. #T2401S010/S100/S500)

PrimeCap is a trademark of Takara Bio Inc.

## Note

This product is not intended for humans or animals *in vivo* applications. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

# Recombinant RNase Inhibitor GMP grade

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