Tris-EDTA Buffer (TE) 10X Powder, pH 7.4

Code No. T9111 Size: 10 pouches

(for 10 L)

Description:

This product is a powder for preparing Tris-EDTA Buffer (TE) commonly used by storing and diluting DNA and RNA. The buffer can easily be prepared by dissolving it with $\rm H_2O.\,$

One pouch is for preparing 1,000 ml of 10X concentrated TE Buffer (pH7.4).

Storage:

Room temperature at desiccated condition.

Specifications:

Format: Powder

Concentration (10X): 100 mM Tris-HCI 10 mM EDTA

Volume: For 1,000 ml pH: 7.4 ± 0.10 at 25° C

Preparation for use:

Transfer one pouch of Tris-EDTA Buffer (TE) 10X Powder, pH 7.4 in a flask or beaker. Add 300 ml of distilled water or deionized water and stir the solution for a few minutes. Adjust the volume to 1,000 ml, and stir until complete dissolution.

Applications:

Storing and diluting DNA and RNA

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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