Takara Stable Competent Cells

Code No. 9132 **Size: 100 μ l x 10**

Supplied Reagents:

pUC19 DNA (0.1 ng/ μ l) **10 μ l**
SOC Medium **1 ml x 10**

Description:

Takara Stable Competent Cells are designed to improve the stability of DNA by disrupting the genes related to homologous recombination. Therefore, they are ideal for cloning of unstable DNA such as the sequences containing poly(A) or repeats. In addition, Takara Stable Competent Cells lack a *fhuA* gene to provide resistance to T1 bacteriophage infection. For cloning, the use of In-Fusion® Snap Assembly Master Mix (Cat. #638943/638944/638947 - 638949) is highly recommended. The competent cells can also be used for blue/white screening using X-Gal* (i.e., β -galactosidase α -complementation) when transformed with pUC-based plasmids.

* X-Gal: 5-Bromo-4-Chloro-3-Indolyl- β -D-Galactoside

Storage: -80°C

Note: Inadequate temperature control may reduce the transformation efficiency. In such cases, check the transformation efficiency using the provided pUC19 DNA. Do not store the competent cells in liquid nitrogen.

SOC Medium:

2%	Tryptone
0.5%	Yeast extract
10 mM	NaCl
2.5 mM	KCl
10 mM	MgSO ₄
10 mM	MgCl ₂
20 mM	Glucose

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Usage:

Cloning of unstable DNA such as the sequences containing poly(A) or repeats

Genotype:

F⁻, *endA1*, *supE44*, *thi-1*, *recA1*, *relA1*, *gyrA96*, *phoA*, Δ *sbcd*, Φ 80d *lacZ* Δ *M15*, Δ *fhuA*, Δ (*lacZYA-argF*) *U169*, Δ (*mrr-hsdRMS-mcrBC*), Δ *mcrA*, λ ⁻

Cell density:

1 - 2 x 10⁹ bacteria/ml

Protocol (When transforming with a plasmid vector):

1. Thaw Takara Stable Competent Cells in an ice bath just before use.
2. After thawing, mix gently to ensure even distribution, and then move 100 μ l of competent cells into a 14 ml round-bottom tube (falcon round tube). Do not vortex.
3. Add no more than 10 ng of DNA for transformation.
4. Place tubes on ice for 30 minutes.
5. Heat shock the cells for exactly 45 seconds at 42°C.
6. Place tubes on ice for 1 - 2 minutes.
7. Add SOC Medium to bring the final volume to 1 ml. SOC Medium should be warmed to 37°C before using.
8. Incubate by shaking (160 - 225 rpm) for 1 hour at 37°C.
9. Plate an appropriate amount of culture on selective medium.*1,2
10. Incubate overnight at 37°C.

*1 If necessary, dilute with the SOC Medium added in step 7.

*2 An appropriate amount should be less than 100 μ l for a 9 cm diameter plate.

Precautions for Use:

- Use no more competent cells than necessary. Transport cells on dry ice/ethanol.
- You may use 1.5 ml microcentrifuge tubes instead of 14 ml round-bottom tubes for transformation (see Protocol Step 2), but this may reduce efficiency.
- When transforming 100 μ l of competent cells, do not use more than 10 ng of purified sample DNA. If you use more than 10 ng of DNA, transformation efficiency may decrease.
- If you change the amount of competent cells, or types of tubes used, it might be necessary to reevaluate the optimal conditions. For example, when using 1.5 ml microcentrifuge tubes, heat shock for 60 seconds at 42°C (see Protocol Step 5).
- Add an appropriate antibiotic to the medium depending on the vector used for cultivation.
- When adding X-Gal to medium, add 20 mg/ml X-Gal (dissolved in dimethylformamide) into 200 - 300 μ l/100 ml agar medium.
- Do not refreeze competent cells once thawed.

References:

- 1) Hanahan, D. *J Mol Biol.* (1983) **166**: 557-580.
- 2) Yanisch-Perron, C., et al. *Gene.* (1985) **33**: 103-119.

Related Products:

In-Fusion® Snap Assembly Master Mix

(Cat. #638943/638944/638947 - 638949)

X-Gal (5-Bromo-4-Chloro-3-Indolyl- β -D-Galactoside) (Cat. #9031)

In-Fusion is a registered trademark of Takara Bio USA, Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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