Vector Map



pLVX-CherryPicker1

Catalog No(s). 632580 (Not sold separately). Sold as part of 632574. Amount

20 µl pLVX-CherryPicker1 Vector (500 ng/µl)

Lot Number Specified on product label.

pLVX-CherryPicker1 Vector Information





Figure 1. pLVX-CherryPicker1 vector map.

 Smal
 Notl

 Apal
 BamHI
 PmeI

 2201
 ACCGCGGGCC
 CGGGATCCGC
 GGCCGCGTTT

 AAACGATGAT
 TGGCGCCCGG
 GCCCTAGGCG
 CCGGCGCAAA



Location of Features

- 5' LTR (5' long terminal repeat): 1–635
- PBS (primer binding site): 636–653
- Ψ (packaging signal): 685–822
- RRE (Rev-response element): 1303–1536
- cPPT/CTS (central polypurine tract/central termination sequence): 2028–2151
- MCS (multiple cloning site): 2206–2235
- CherryPicker (mCherry-transferrin receptor membrane-anchor domain fusion): 2236–3162
- P_{PGK} (phosphoglycerate kinase promoter): 3176–3684
- Puro^r (puromycin resistance gene): 3705–4304
- WPRE (woodchuck hepatitis virus posttranscriptional regulatory element): 4318–4909
- 3' LTR (3' long terminal repeat): 5112–5748
- pUC origin of replication: 6217–6890 (complementary)
- Amp^r (ampicillin resistance gene; β-lactamase): 7035–8031 (complementary)

NOTICE TO PURCHASER:

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements listed below and described on the product's web page at http://www.clontech.com. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.